INHEALTH[®]

TECHNOLOGIES

Blom-Singer[®] Indwelling Voice Protheses

for Primary, Secondary,
Replacement TEP Procedures

MEDICAL PROFESSIONAL Instructions For Use

Blom-Singer[®]

voice restoration systems

R1

37640-01J

37640-01J Effective June 2018 / Gültig ab Juni 2018 / Efectivo a partir de junio de 2018 / Date d'entrée en vigueur : juin 2018 / Valido da giugno 2018 / Met ingang van juni 2018 / Vigente em junho de 2018 / V platnosti od června 2018 / Wchodzi w życie w czerwcu 2018 r.

Blom-Singer and **InHealth Technologies** are registered trademarks of **Freudenberg Medical, LLC**

Patent(s): www.inhealthpatents.com

TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Índice / Obsah / Spis treści

ENG	Indwelling Voice Protheses	4
DEU	Verweil-Stimmprothesen	18
SPA	Prótesis fonatorias fijas	34
FRA	Prothèses vocales à demeure	50
ITA	Protesi fonatorie fisse	66
NLD	Verblijfstemprothesen	82
POR	Próteses vocais fixas	98
CES	Vnitřní hlasové protézy	114
POL	Stałe protezy głosowe	128

Blom-Singer®

voice restoration systems

MEDICAL PROFESSIONAL Instructions For Use

Gebrauchsanweisung für MEDIZINISCHE FACHKRÄFTE
Instrucciones de uso para PROFESIONALES MÉDICOS
Instructions d'utilisation destinées aux PROFESSIONNELS DE SANTÉ
Istruzioni per l'uso destinate al PERSONALE SANITARIO
Gebruiksaanwijzingen voor de MEDISCHE PROFESSIONAL
Instruções de uso para PROFISSIONAIS MÉDICOS
Návod k použití pro ZDRAVOTNICKÝ PERSONÁL
Instrukcja użytkowania dla PERSONELU MEDYCZNEGO

BLOM-SINGER® INDWELLING VOICE PROSTHESES

for Primary, Secondary, Replacement TEP Procedures

This Blom-Singer medical device may be used only once by one patient. It may not be reused. See Warnings, Precautions and Complications.

Do not use the device if it becomes contaminated; discard and replace it with a newly packaged prosthesis.

INDICATIONS

The Blom-Singer Indwelling Voice Prosthesis is indicated for tracheoesophageal voice restoration following total laryngectomy, when placement, or replacement, of an indwelling voice prosthesis is performed by a qualified, trained medical professional.

CONTRAINDICATIONS

The Blom-Singer indwelling voice prosthesis is a medical product that should be used only by a qualified, medical professional with experience and training in its use and care.

HOW SUPPLIED

The Blom-Singer Indwelling Voice Prosthesis is supplied in sterile or non-sterile options, as a complete voice restoration system for primary, secondary, or replacement procedures. Sterile voice prostheses are sterilized by ethylene oxide. Included in each package is a Voice Prosthesis Cleaning System for the patient, two (2) gel caps, flange introducer and inserter.

PRODUCT DESCRIPTION

Please refer to the diagrams located at the back of this instruction manual.

Nonsterile product is not intended for primary placement procedures.

The Blom-Singer voice prosthesis is made of silicone, a

biocompatible material used in many medical products.

The gel cap used for insertion of the prosthesis is made of a vegetable base gelatin. Individual patient reactions to these materials may occur.

The Blom-Singer voice prosthesis, (diagram 1), consists of a one-way silicone flap valve (diagram 1a), an esophageal retention flange (diagram 1b), a body that holds the valve assembly (diagram 1c), a tracheal retention flange (diagram 1d), and a neck strap (diagram 1e).

WARNINGS AND PRECAUTIONS

The Blom-Singer indwelling voice prosthesis is a medically prescribed device. Initial sizing and instructions for use must be provided by a qualified medical professional trained in the use of this particular device.

The qualified, trained medical professional is to instruct patients on the use and care of this device and to provide patients with the Blom-Singer instructions for use supplied with the product.

If the prosthesis is not functioning properly, the patient should have it evaluated by a qualified, trained medical professional as soon as possible.

If there are tears, cracks, or structural damage to the prosthesis, cleaning device or accessories, discontinue use.

Handle the device and the gel cap loading tool carefully to avoid damaging the prosthesis.

When using a sterile voice prosthesis, aseptic handling technique must be used. All indwelling voice prostheses must be carefully removed from their packaging and handled in a manner to prevent contamination of the devices.

Use only the specially sized Blom-Singer Gel Caps.

Do not use petroleum-based products such as Vaseline® to lubricate the voice prosthesis as these products can damage the device.

Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. It is

recommended that the voice prosthesis and TEP be evaluated by the clinician at regular intervals.

Voice Production

To prevent post-operative complications, the patient should not begin speaking with the voice prosthesis until the clinician has indicated that it is safe to do so.

The lumen of the voice prosthesis must be kept clear of blockage in order for it to function properly to allow the patient to voice.

In some users the inability to relax the muscles of the throat may account for their inability to speak fluently and, to speak with minimal effort. This problem requires professional assessment.

Patients requiring postoperative radiation may have transient interruption of voice in the third or fourth week of treatment. The prosthesis can be left in place or a catheter can be inserted as a stent until radiation treatments are completed.

Voice Prosthesis Dislodgement

Care must be exercised during voice prosthesis insertion or removal to avoid injury to the TEP or accidental displacement of the voice prosthesis which could result in aspiration of the voice prosthesis into the trachea. Should aspiration occur, the patient should attempt to cough the voice prosthesis out of the trachea. Further medical attention may be necessary if coughing the voice prosthesis out is unsuccessful.

Confirm gel cap dissolution and deployment of the esophageal retention flange per procedures provided below to ensure device is securely retained in TEP.

If the voice prosthesis is accidentally dislodged from the puncture, the patient should immediately place a Blom-Singer Puncture Dilator or suitable device, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. The voice prosthesis should be reinserted within 24 hours.

Foreign objects, other than the supplied cleaning devices, should never be inserted into the voice prosthesis by the patient.

Inserting objects other than the supplied cleaning devices may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or aspirated.

Care should be exercised when inserting or removing cleaning devices to avoid accidental displacement of the voice prosthesis which could result in aspiration of the voice prosthesis.

Voice Prosthesis Leakage

When the flap valve fails to close completely a few drops of fluid may pass back through the valve from the esophagus to the trachea, which may cause coughing. Recurrent leakage of the voice prosthesis should be evaluated.

Gentle handling and pressure should always be used when cleaning the voice prosthesis to avoid deterioration of the device and damaging internal components which could cause leakage.

Microbial Growth

Most users have microbial growth in their mouth, throat, or esophagus; this is a commonly occurring condition. Microbial growth deposits can cause valve deformation and failure, i.e., liquid leakage through the voice prosthesis and an increase in the pressure necessary to voice. Replacement of the voice prosthesis may be required.

Device Removal

The attached neck strap on an indwelling voice prosthesis should only be removed after the esophageal retention flange placement within the esophagus has been verified. Never attempt to insert or reinsert an indwelling voice prosthesis that has the neck strap removed.

Never remove one indwelling voice prosthesis and insert another voice prosthesis without first dilating the TE puncture and re-measuring the tract to confirm correct voice prosthesis length.

The patient should never attempt to remove the device or allow anyone other than a qualified, trained medical professional to remove the indwelling voice prosthesis.

The indwelling voice prosthesis is not a permanent device and requires replacement periodically.

The indwelling voice prosthesis may be left in place in the TE puncture until it ceases to function correctly, that is, until it has persistent leakage or is not providing adequate voice for speech.

Many clinicians recommend routine follow-up evaluations not longer than 6 months from initial placement.

Removal of the indwelling voice prosthesis should only be done by grasping the outer tracheal flange of the device securely with a locking hemostat. Pull gently and firmly until the prosthesis is fully withdrawn.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone prostheses of the Blom-Singer type. They include: stoma contamination or sepsis, which may require removal of the voice prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; occasional extrusion of the prosthesis, requiring replacement after dilation of the TE puncture and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the voice prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TE puncture; intractable leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia, tearing or other damage to the voice prosthesis from improper use; microbial growth deposits causing leakage through the voice prosthesis or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus.

INSTRUCTIONS FOR USE

The following instructions are provided for primary and secondary/replacement procedures.

Always measure the TE puncture tract **before** inserting a voice prosthesis to confirm correct voice prosthesis length.

Primary Voice Prosthesis Placement

The following procedure is provided by the late Ronald C. Hamaker, M.D., F.A.C.S., as recommended instructions for primary placement of the voice prosthesis at the time of the tracheoesophageal puncture during laryngectomy.

Prior to use of the voice prosthesis, check the valve mechanism to be sure it is intact and working properly. The flap valve, (diagram 8a), should close flat against the seating surface inside the voice prosthesis.

Please refer to the diagrams located at the back of this instruction manual.

1. Grasp the neck strap of the sterile indwelling voice prosthesis with a right angle hemostat (diagram 8b).
2. Insert the tip of the hemostat with neck strap into the esophagus. Push the tip against the anterior wall 1.5 cm. from the cut edge of the trachea. Cut down horizontally on the tip through the trachea (diagram 9).
3. The hemostat tip with the neck strap is gently pushed through the tracheoesophageal wall into the trachea (diagram 10).
4. The neck strap is grasped with a second hemostat and the voice prosthesis is pulled through the puncture site (diagram 11).
5. The voice prosthesis is positioned in the tracheoesophageal puncture (diagram 12).
6. The neck strap on the voice prosthesis is detached and the wound is closed. Voicing is deferred until two weeks postoperatively to avoid the possibility of an air pressure induced fistula (diagram 13).

Secondary Voice Prosthesis Placement

The following procedure is provided by Eric D. Blom, Ph.D. as recommended instructions for secondary voice prosthesis placement, or replacement of the voice prosthesis, using the Blom-Singer® Gel Cap insertion method.

Please refer to the diagrams located at the back of this instruction manual.

Note: All secondary procedures must only be done using a bright light to illuminate the tracheostoma and puncture. Insertion or removal of the voice prosthesis in the TE puncture should only be done with a bright light focused directly on the stoma.

Use of gloves and protective eye wear is recommended for the following procedure.

Have the following Blom-Singer items prepared for insertion of the voice prosthesis:

- Blom-Singer indwelling voice prosthesis
- Blom-Singer Gel Caps
- Blom-Singer Dilator-Sizer or Blom-Singer Voice Prosthesis Sizer, and Blom-Singer Puncture Dilator
- Blom-Singer voice prosthesis inserter
- Flange introducer

Dilate and Measure

Please refer to the Blom-Singer Dilation/Sizing System Instructions For Use for complete product details.

Optional: The Blom-Singer Sizer or Puncture Dilator may be used as an alternate procedure. Refer to the instructions for use supplied with the product.

1. Dilate the TEP with a Blom-Singer Dilator-Sizer of the appropriate diameter. Gently insert the Dilator-Sizer into the puncture completely to the neck strap (diagram 2) and leave in place for five minutes.
2. Carefully withdraw the Dilator-Sizer until very gentle resistance is detected indicating that the esophageal wings are seated against the anterior wall of the esophagus (diagram 2a). Note: Pulling too firmly on the Dilator-Sizer may squeeze or compress the TE tract resulting in an erroneous measurement that indicates it is shorter than it actually is.
3. Note the marker that is closest to the entrance of the TE puncture. Read the size number on the Dilator-Sizer and select a Blom-Singer Indwelling Voice Prosthesis of corresponding size. If the measurement is between two markers, select the longer one.
4. Fully insert the Dilator-Sizer and fasten its neck strap with medical adhesive tape on to the patient's peristomal skin. Leave the Dilator-Sizer in the puncture to dilate 5 minutes, or longer if necessary, in order to allow for easy insertion of the selected indwelling voice prosthesis. Prior to insertion of the Blom-Singer Indwelling Voice Prosthesis, remove the tape

from the Dilator-Sizer neck strap, grasp the dilator securely, and carefully withdraw it from the puncture. Immediately insert the indwelling voice prosthesis (see below).

Voice Prosthesis Preparation

Note: Prior to use of the indwelling voice prosthesis, check the valve mechanism to be sure it is intact and working properly. The flap valve, (diagram 1a), should close flat against the seating surface inside the voice prosthesis.

Prepare the voice prosthesis for insertion by utilizing the Blom-Singer Gel Cap insertion method. Use only the specially sized gel caps provided with this product.

The following steps describe the hand folding method of gel cap insertion.

1. Remove the voice prosthesis from the package. Always handle the voice prosthesis with clean, gloved hands. Be certain that your hands and the voice prosthesis are completely dry before loading the gel cap. Failure to do so may cause premature dissolution of the gel cap.
2. Remove the gel cap insertion components from the package: gel cap, flange introducer, and inserter tool.
3. Fold and load the esophageal retention flange of the voice prosthesis into the gel cap, (shorter end of the capsule) according to the following steps. Discard the gel body, (longer end of capsule may be printed), see figure below.

3.1. Fold the esophageal retention flange tightly in half toward the center of the device (diagram 4a).

3.2. Fold the two outer edges of the esophageal retention flange tightly against the folded flange (diagram 4b) and hold the prosthesis in this folded position.

3.3. Insert the folded portion of the esophageal retention flange into the gel cap (diagram 4c). Gently tuck and push the folded esophageal retention flange as far as possible into the gel cap.

3.4. Use the flange introducer to gently tuck the

remaining portion of the folded esophageal retention flange completely into the gel cap (diagram 4d). Do not use any sharp or serrated instruments to avoid damaging the prosthesis.

3.5. The esophageal retention flange should be folded in a forward position and completely enclosed inside the gel cap (diagram 4e).

Prosthesis Insertion

Warning: Caution should be exercised to avoid accidental displacement of the voice prosthesis which could result in aspiration of the voice prosthesis into the trachea. Should aspiration occur, the patient should attempt to cough the voice prosthesis out of the trachea. Further medical attention may be necessary if coughing the voice prosthesis out is unsuccessful.

1. Place the voice prosthesis on the inserter and lock the neck strap onto the safety peg (diagram 3).
2. Caution the patient not to swallow as you remove the previously inserted dilator-sizer to avoid aspiration of saliva into the trachea prior to inserting the voice prosthesis.
3. Apply a light coating of water-soluble lubricant (oil-free) to the tip of the gel-capped end of the voice prosthesis and immediately insert the voice prosthesis completely into the TE puncture by aligning the tip of the voice prosthesis partially in the puncture with the neck strap oriented upwards (diagram 5). **Note:** Do not use petroleum-based products such as Vaseline® to lubricate the voice prosthesis as these products can damage silicone.
4. Hold the voice prosthesis in this position of complete insertion for at least three minutes. This allows ample time for the gel cap to dissolve and release the esophageal retention flange (diagram 1b) within the esophagus. Patients should be instructed to swallow their saliva or warm water to facilitate dissolving the gel cap.

Esophageal Flange Deployment of Dual Valve Indwelling Voice Prosthesis

To accommodate esophageal flange deployment of the Dual Valve Indwelling Voice Prosthesis, the inserter is designed so that it does not fully insert into the voice prosthesis (diagram 7). Within three minutes the inserter should visibly move forward into the voice prosthesis; the clinician should see this forward movement indicating that the gel cap has dissolved and that the esophageal retention flange has deployed.

Caution: If the voice prosthesis does not insert easily on the first attempt, do not continue to try to insert it. Instead, re-insert the puncture dilator, of the appropriate diameter, for a few minutes to dilate the TE tract, and then try again.

5. Detach the voice prosthesis neck strap (diagram 1e) from the safety peg of the inserter. Place a finger against the neck strap and carefully withdraw the inserter from the voice prosthesis with a twisting motion. **Warning:** Caution should be exercised when removing the inserter to avoid accidental displacement of the voice prosthesis which could result in aspiration of the voice prosthesis. Secure the voice prosthesis with medical adhesive tape across the neck strap. It is recommended to orient the prosthesis with the neck strap up.
6. Have the patient take a breath, occlude the stoma, and produce voice. Ensure there is no leakage through the voice prosthesis by having the patient drink water.

Confirmation of Esophageal Retention Flange Deployment

The following procedure for confirming deployment of the esophageal retention flange is recommended by Eric D. Blom, Ph.D.

After insertion of the voice prosthesis in the puncture tract and completion of the 3-minute waiting period for gel cap dissolution and esophageal retention flange deployment, repeatedly rotate the prosthesis on the inserter while it is in the puncture. If the esophageal retention flange has deployed and has properly seated against the mucosa of the anterior esophageal wall, the voice prosthesis on the inserter will rotate easily in the puncture

tract. If the gel cap has dissolved and incorrectly deployed the esophageal retention flange inside the puncture tract rather than inside the lumen of the esophagus, considerable resistance to voice prosthesis rotation will be felt.

A second method for confirming gel cap dissolution and esophageal retention flange deployment is to assess the capacity to produce tracheoesophageal voice. If the voice prosthesis is correctly positioned and the esophageal retention flange deployed, airflow diverted from the trachea by digital occlusion of the tracheostoma will open the flap valve and pass up the esophagus to initiate fluent voice.

Proper Valve Function Confirmation of the Dual Valve Indwelling Voice Prosthesis

Use a small wax curette or equivalent gauge instrument to gently open the tracheal valve while the patient swallows their saliva or water. The clinician should see the properly positioned esophageal valve at the other end and no leakage.

If the proper deployment and positioning of the esophageal retention flange cannot be confirmed by one or both of the methods described above, the clinician may need to determine gel cap dissolution and esophageal retention flange deployment by manually opening, with gentle contact, the one-way flap valve in the esophageal end of the voice prosthesis.

Note: The following method is not applicable to the Dual Valve Indwelling Voice Prosthesis.

Direct a bright light into the prosthesis to visualize the flap valve at the proximal end. **Slowly** and **carefully** advance the shaft end (non-cotton tip) of a long cotton swab into the voice prosthesis until you observe that **gentle contact** opens the flap valve. A flap valve that does not open with gentle contact indicates that the undeployed retention flange or an inverted flap valve is preventing valve opening.

If the assessment methods described above fail to provide satisfactory verification of correct esophageal retention flange deployment, remove the voice prosthesis, dilate and re-measure the puncture tract, and repeat the voice prosthesis insertion and confirmation process.

Neck Strap Detachment

The indwelling voice prosthesis is designed to include the optional detachment of the neck strap following confirmation that the voice prosthesis esophageal retention flange is **securely** positioned within the esophageal lumen against the anterior wall of the esophagus. Removal of the neck strap is an **optional** procedure.

Once deployment of the esophageal retention flange has been confirmed, detach the neck strap from the tracheal retention flange at the area of reduced neck strap width, where it meets the tracheal flange, by cutting carefully with scissors (diagram 6).

The neck strap should not be detached from the indwelling voice prosthesis if:

- use of a laryngectomy tube may potentially dislodge the prosthesis;
- a patient has previously been wearing a voice prosthesis that is larger in diameter than the indwelling voice prosthesis currently being inserted;
- the diameter of the TE puncture is dilated significantly larger than the voice prosthesis diameter; or,
- the health of the stoma tissue or physical condition of the patient is in question.

Note: It is important **not** to detach the neck strap from the voice prosthesis in the instances cited above. Tape the neck strap to the peristomal skin and observe security of voice prosthesis retention for an extended period of time.

Removal of the Prosthesis

The indwelling voice prosthesis is not a permanent device and requires replacement periodically. The indwelling voice prosthesis may be left in place in the TE puncture until it ceases to function correctly, that is, until it leaks or is not providing adequate voice for speech. Many clinicians recommend routine follow-up evaluations not longer than 6 months from initial placement.

Removal of the indwelling voice prosthesis should only be

done by grasping the outer tracheal flange of the device securely with a locking hemostat. Pull gently and firmly until the voice prosthesis is fully withdrawn.

Insert a dilator-sizer of the appropriate diameter, and tape it in position for five minutes prior to inserting a new indwelling voice prosthesis.

Warning: The patient should **never** attempt to remove the device or allow anyone other than a qualified, trained medical professional to remove the indwelling voice prosthesis.

Never remove one indwelling voice prosthesis and insert another voice prosthesis without first dilating the TE puncture and remeasuring the tract to confirm correct voice prosthesis length.

VOICE PROSTHESIS CLEANING AND CARE INSTRUCTIONS

The following instructions should be made clear to the patient as part of the routine care and cleaning of the indwelling voice prosthesis.

The following information is a brief overview of the proper care and maintenance of the voice prosthesis. Please refer to the Voice Prosthesis Cleaning System instructions for use supplied with this product for complete instructions on how to use the cleaning devices.

The purpose of cleaning the voice prosthesis is to remove blockage of its lumen that impairs the patient's ability to voice.

The cleaning device provides a means for careful cleaning of the voice prosthesis while in the patient's TE puncture.

Caution: Use **only** supplied cleaning device. **Do not** insert foreign objects into voice prosthesis. Inserting objects other than the supplied cleaning device may cause dislodgment of the voice prosthesis or its components.

The user's hands should always be washed with soap and water prior to using the cleaning system.

Cleaning of the indwelling voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the prosthesis is clearly visible.

Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris which can be aspirated into the airway.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly.

SPECIAL ORDER PRODUCTS

If this instruction manual accompanies a Special Order Product, there may be differences in the physical characteristics between the product and the product descriptions described in this instruction manual. These differences will not affect the safety or efficacy of the special order product. Special Order products are nonreturnable.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday — Friday, 9:30 am — 7:00 pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service.

International

Please contact our customer service department for a distributor referral.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

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












Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.

-
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
 11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
 12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
 13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
 14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
 15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
 16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

INHEALTH®

TECHNOLOGIES

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

CE 0344

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands