FAHL

NEBUJUNIOR®

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CONTENT

EN

1.0	Introduction	25	7.0	Maintenance/reprocessing for reuse in another patient	36
1.1	Notes on the instructions for use	25	7.1	Basic instructions	36
1.2	Function	26	7.2	Reprocessing for reuse	36
1.3	Intended use	26	7.2.1	Air filter	36
1.4	Scope of delivery	27	7.2.2	Nebulizer	36
1.5	Transport and storage	27	7.2.3	Air tube	36
1.6	Legend of the Pictograms	27	7.2.4	Attachments and Adapters	36
			7.2.5	Device surface	36
2.0	Safety information	28	8.0	Troubleshooting	37
3.0	Set-up and start-up	30	9.0	Accessories, consumables	38
3.1	Operating elements	30	10.0	Technical data	39
3.2	Start-up	30	11.0	Disposal	40
4.0	General operation	31	12.0	Instructions with regard to EMC	41
4.1	Connecting the nebulizer	31	12.1	Guidelines and manufacturer's declaration – electromagnetic immunity	41
4.2	Connecting the air tube	31	12.2	Guidelines and manufacturer's declaration – electromagnetic immunity for the NEBUJUNIOR® Penguin	42
4.3	Connecting attachments	31	12.3	Guidelines and manufacturer's declaration - electromagnetic emissions	43
5.0	Operation	32	12.4	Recommended safety distances between portable and mobile HF telecommunication devices and the NEBUJUNIOR® Penguin	43
5.1	Nebulizer	32	13.0	Legal notices	44
5.2	Attachments	32			
6.0	Cleaning, disinfection, reprocessing	33			
6.1	Basic instructions	33			
6.1.1	Air filter	33			
6.1.2	Air tube	34			
6.1.3	Nebulizer and attachments	34			
6.1.4	Device surface	34			
6.2	Disinfection	35			
6.3	Recommended instrument disinfectants	35			
6.4	Recommended surface disinfectants	35			

To request further information, accessories, consumables or spare parts, contact:

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1.0 INTRODUCTION

1.1 NOTES ON THE INSTRUCTIONS FOR USE

These instructions for use contain important information as to how you can operate the NEBUJUNIOR® Penguin safely, properly and effectively. It helps avoid hazards and reduce repair costs and outage times, thus enhancing, amongst other features, the reliability and prolonging the service life of the device. They not only serve for briefing persons who operate the device, but are also intended to be used as reference.

The instructions for use must always be kept available close to the device.

Technically correct application ensures the operational safety and the functionality of the NEBUJUNIOR® Penguin and is therefore indispensable, besides cleaning performed on a regular basis.

Repairs and professional reprocessing for reuse may only be performed by the manufacturer or by an authorised specialist company. By using original spare parts, you have the warranty that the operational safety, availability for use as well as the value of your NEBUJUNIOR® Penguin will be preserved.

- Read Section 2.0 „Safety information“ before taking the device into operation for the first time in order to avoid possible hazardous situations.
- These instructions for use correspond to the current state of the technical safety standards at the time of going to print.
- Reprinting – including in excerpts – is subject to the written consent of Andreas Fahl Medizintechnik Vertrieb GmbH.
- Subject to alterations. Errors excepted.

These instructions for use are valid for the following device:

NEBUJUNIOR® Penguin (REF 50400)

In connection with the accessory parts specified in Section 9.0, as well as for all devices of identical construction with special voltage.

Please keep this document for later use!

1.0 INTRODUCTION

EN


Figure 1


Figure 2

NEBUJUNIOR® Penguin

- ❶ Air compressor
- ❷ Ventilation slits
- ❸ Air filter compartment/ Air filter insert
- ❹ Air connection
- ❺ Mains connection cable
- ❻ On/Off button

1.2 FUNCTION

- The NEBUJUNIOR® Penguin is a modern Inhaler device for compressed-air operated humidification of breathing air.
- All medications approved for use in inhalation therapy, as well as other inhalants, may be nebulized using this device. Your physician will be happy to assist you in selecting medications appropriate to your needs.
- Using the supplied attachments, inhalation is alternatively possible through the mouthpiece, the nosepiece or the inhalation mask.
- Special accessories include the LARYVOX® mask, a tracheostoma mask for tracheostomy patients.
- A PEP cap is also available.
- The NEBUJUNIOR® Penguin can easily be disassembled into the individual components for maintenance and cleaning.

1.3 INTENDED USE

Name:
NEBUJUNIOR® Penguin

Main function:
Humidification of breathing air through inhalation of steam and medications and other inhalants approved for use in inhalation therapy.

Medical indication / application:
Humidification of the upper and lower airways.

Specification of the main functions:
Prevention and treatment of diseases of the upper and lower airways. The aerosol is produced using compressed air. This enables humidification of the airways.

Application organ:
Upper and lower airways (oral cavity, nasopharyngeal space and bronchial system).

Duration of application:
Temporary patient use for the purposes of inhalation.

Application environment:
The application environment includes hospitals, doctors' practices, nursing facilities and use at home. Where medication is also used, a physician must determine the respective type and dosage. Application must only be carried out by medically trained and briefed persons. Children and patients who require special assistance should only carry out inhalation under professional guidance.

Contraindication:
The NEBUJUNIOR® Penguin must not be used:

- if it is damaged,
- by unsupervised children,
- in damp rooms or in the shower/bath,
- outdoors.

The NEBUJUNIOR® Penguin is an active medical device.

Reprocessing:
The product is intended for repeated use. The accessories are partly reusable. For information on reprocessing, cleaning and disinfection, refer to instructions for use (see Section 6.0).

The NEBUJUNIOR® Penguin must **not** be used:



- if it is damaged, or in combination with medications and devices that have not been approved for use in inhalation therapy;
- in damp rooms, in the shower/bath, or outdoors, or by unsupervised children.

1.0 INTRODUCTION

1.4 SCOPE OF DELIVERY

- This device was subjected to a thorough function check and was carefully packaged before being shipped. Please, nevertheless check the content of the packaging for completeness by comparing it with the scope of delivery immediately upon receipt. Besides the basic unit, the scope of delivery comprises the following:

NEBUJUNIOR® Penguin

	Nebulizer
	Air tube
	Air filters
	Inhalation mask for adults
	Inhalation mask for children
	Nosepiece
	Mouthpiece
	Adapter for 15 mm-connector (UNI)
	PEP cap
	Carrying bag
SPECIAL ACCESSORIES (not included in the scope of delivery)	
	LARYVOX® mask Tracheostoma mask for tracheotomised patients

1.5 TRANSPORT AND STORAGE

- The device may only be transported in the enclosed transport bag.
- Please promptly document and report transport damages.
- Following transport at temperatures below freezing point, the device must be left to stand at room temperature for up to six hours before start-up. If the device is not acclimatised, it must **not** be operated.

- Ambient conditions:**

Transport / storage:	-20 to +70°C; 10 to 95 % air humidity without condensation
Operation:	+10 to +40°C; 30 to 75 % air humidity without condensation

1.6 LEGEND OF THE PICTOGRAMS

Observe instructions for use!

Device of protection class II

Type BF degree of protection

Serial number

Catalogue number

Manufacturer

Please observe the laws and regulations which apply to the disposal of electrical and electronic devices in the respective countries.

Attention, important information

IP classification

2.0 SAFETY INFORMATION

EN

GENERAL SAFETY INSTRUCTIONS

Please read these instructions for use through carefully before starting up the NEBUJUNIOR® Penguin.

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- no original accessory parts and spare parts are used,
- the instructions regarding proper use in these instructions for use are disregarded,
- assembly, readjustments, modifications, extensions and repairs have not been carried out by persons and companies who have been authorised by Andreas Fahl Medizintechnik-Vertrieb GmbH.

The NEBUJUNIOR® Penguin conforms to the accepted rules of technology and the provisions of the Medical Devices Act.

No warranty claims are accepted in the event of damages or malfunctions that have arisen due to the use of accessories or consumables from external companies.

Only transparent tubes must be used.

The device must not be operated in environments exposed to explosion hazard and enriched with oxygen.

RISK OF INJURY!

The NEBUJUNIOR® Penguin is designed for humidification of the upper and lower airways. Only medications or other inhalants approved for use in inhalation therapy may be used.

The device must not be operated in areas exposed to splash water or in zones exposed to explosion hazard.

The NEBUJUNIOR® Penguin may only be used by persons who have been briefed in the medical application.

To disconnect from the mains power supply, pull the mains plug out of the wall socket:

- before cleaning and performing maintenance on the device,
- each time before leaving the room.

Do not pull the plug out of the wall socket by pulling on the cable!

Never touch the plug or the cable with wet hands.

CAUTION!

There is a risk of strangulation due to the long air tube or the power supply cable.

EQUIPMENT DAMAGE!

The device must not be taken into operation:

- **if the power supply cable or the connection socket are damaged,**
- **if the device has fallen down,**
- **if the device shows obvious safety defects.**

Clean the device and return it for repair via the dealer from whom you have purchased the NEBUJUNIOR® Penguin.

The device must be checked at regular intervals for proper function and technical safety defects, e.g. plug contacts, housing, etc.

The ambient conditions specified in Section 1.5 ,Transport and storage' must be observed.

Never immerse the device in water, even when the device is in non-activated state.

Before starting up the device, the device itself and accessories must be checked for damages. Ventilation slits should not be obstructed. The device must be checked for proper function before use.

No liquid must be allowed to penetrate into the device. If liquid has penetrated into the device, then the device must only be taken into operation again after it has been checked by the customer service. Damaged accessories must be replaced immediately.

The device must be set up on a firm, level surface. If the device is set up on an uneven surface (e.g. mattress, cushion, upholstered chair), there is a risk that the device may tip.

Before the device is connected, it must be checked whether the mains voltage and frequency indicated on the device are the same as the voltage and frequency of the supply network.

Only use proper mains connections and extension cables. Avoid moisture on plug and switch unit.

3.0 SET-UP AND START-UP

3.1 OPERATING ELEMENTS

- ⑥ On/Off button
- ⑦ Mains connection cable


Figure 3

3.2 START-UP

- The NEBUJUNIOR® Penguin is delivered ready for operation.
- Lift the device out of the packaging. Carry out a visual inspection for external damage.
- Check whether the voltage values indicated on the power rating plate agree with the voltage supplied.
- Always place the device on a level, secure surface.
- You must observe the safety instructions in Section 2.0 before starting up the device for the first time.
- Following transport at temperatures below freezing point, the device must be left to stand at room temperature for up to six hours before start-up. If the device is not acclimatised, it must not be operated.
- Only plug the NEBUJUNIOR® Penguin into a 230 V wall socket.
- Do not place the NEBUJUNIOR® Penguin on cloths or blankets.
- Ventilation slits should not be obstructed.

4.0 GENERAL OPERATION

Please read the preceding section before going on to read these instructions for use!

IMPORTANT SAFETY INSTRUCTIONS!

Caution: Do not perform inhalation in the airways until you have been briefed by specialist staff.

The nebulizer should always be held in an upright position (Figures 17-20).

Only use the device and accessories if they are clean and disinfected.

Unless otherwise prescribed by the physician, inhalation should be performed 2 – 3 times daily for approx. 10 – 15 minutes.

Once phlegm has loosened, it should be coughed up, or in the case of tracheostomy patients, suctioned using a suctioning device and a sterile suction catheter.

Draughts and cold temperatures should be avoided immediately after inhalation.

4.1 CONNECTING THE NEBULIZER

- Pull off the top section of the nebulizer by twisting both parts (anti-clockwise) gently against each other (Figure 4).
- In the lower section of the nebulizer (Figure 5), pour in the exact quantity of the medication/inhalant in accordance with the instructions provided by the manufacturer or the attending physician.
- Observe the maximum capacity of 6 ml.
- Place the top section of the nebulizer back onto the lower section. Connect the nebulizer by twisting both parts (clockwise) gently against each other. Ensure that these fit together securely (Figure 6).

Figure 4


Figure 5

Figure 6


Figure 7


Figure 8


Figure 9


Figure 10


Figure 11


Figure 12

4.2 CONNECTING THE AIR TUBE

- Check whether the air tube is bent or kinked. Connect the air tube (Figure 7) with the lower opening of the nebulizer.
- Then connect the air tube with the air connection of the device (Figure 8).

4.3 CONNECTING ATTACHMENTS

- Attach the mouthpiece (Figure 9) to the side connection of the nebulizer. Please also use the adapter (Figure 10) included in the delivery in order to make use of the nosepiece (Figure 11), the inhalation mask (Figure 12) or the LARYVOX® mask (Figure 13).
- In order to increase breathing resistance, you can place the PEP cap included in the delivery on to the expiration valve on the top section of the nebulizer (Figure 14).


Figure 13

Figure 14

5.0 OPERATION

EN

5.1 NEBULIZER

- Hold the nebulizer in an upright position where possible.
- Two applications are possible: continuous inhalation and temporary inhalation. For continuous inhalation, push the locking device (Figure 15 a/b) downwards in the direction indicated by the arrow in order to secure the blue lever for the desired duration of nebulisation. For temporary inhalation, hold down the blue lever.
- To perform inhalation, sit in a relaxed pose and breath in and out evenly.
- To allow your airways to become accustomed to the aerosol, do not breathe in too deeply during the first breaths you inhale.
- In order to increase breathing resistance, you can place the PEP cap on to the expiration valve on the top section of the nebulizer (Figure 16).


Figure 15a


Figure 15b


Figure 16

5.2 ATTACHMENTS

- When using the mouthpiece, enclose it fully with your lips (Figure 17).
- During inhalation, the LARYVOX® mask (special accessories) must be connected to the nebulizer (Figure 12 - 13) using the adapter.
- Press the inhalation mask/tracheostoma mask gently over your face/tracheostomy/nose (Figure 18 - 20), relax, and breath in.
- During inhalation, the nosepiece must be connected to the nebulizer using the adapter. Then hold the nosepiece to the nose (Figure 20).
- When breathing out, the inhalation mask/tracheostoma mask/nosepiece should be held a small distance (the width of a finger) away from your face/tracheostomy/nose.


Figure 17


Figure 18


Figure 19


Figure 20

IMPORTANT SAFETY INSTRUCTIONS FOR THE DEVICE

The device must not be operated without the air filter!

Therefore, have at least one spare air filter ready!

Check each time before use whether the air filter is clean and dry. Moist or contaminated air filters must be replaced by new air filters.

For your own safety, use disposable gloves when replacing the air filter!

The air filter is a single-patient product and must be replaced at least every 4 weeks (depending on frequency of use).

IMPORTANT SAFETY INSTRUCTIONS FOR THE DEVICE

Risk of fatal injury due to electric shock.

Prior to cleaning and disinfection, pull the mains plug out of the wall socket.

Careful handling of the inhalation device has a decisive impact on the reliability and safety of the device. A regular visual inspection of the device for external damage and contamination is indispensable.

The cleaning measures described in the following must be carried out even when the device is not being used regularly. Regular replacement of the air filters for the NEBUJUNIOR® Penguin at 4-week intervals and function checks at intervals of not more than 3 months are mandatory in order to ensure the functionality of the device.

These hygiene measures are necessary measures to protect patients and users and to maintain the functional reliability of the inhalation device.

If the device is to be used for another patient, professional reprocessing at the manufacturer's establishment or at an authorised specialist company in accordance with the Medical Devices Operator Ordinance (MPBetreibV), the Medical Devices Act (MPG) and the Federal Association for Medical Technology (BV-Med) guideline is required.

CAUTION!

The inhalation device is a single-patient product and may only be used for one patient, since otherwise there is a risk of infection due to viruses, bacteria or germs. Infections can cause the spread of life-threatening diseases. Following inhalation, the device and accessories must be cleaned and disinfected. This can minimise the risk of infection. This inhalation device can be reprocessed for reuse by another patient. When doing this, observe the instructions in this regard provided in the instructions for use, particularly in Section 7. (Maintenance / reprocessing for reuse in other patient). The manufacturer accepts no liability for unauthorised and incorrect reprocessing of this device by external companies!

6.1 BASIC INSTRUCTIONS

The inhalation device must be cleaned in accordance with the following cleaning and maintenance instructions.

The following instructions are the result of long years of experience. Shorter cycles may be required, depending on the application and the users' experience.

- A regular visual inspection of the device for external damage and contamination is indispensable.
- Depending on the use of the device, e.g. for nursing purposes, additional protective measures must be taken before beginning inhalation, to ensure that the safety of the patient, the user and third parties is not jeopardised. For example, hands should be thoroughly disinfected.
- Disconnect the power supply cable before cleaning.
- The cleaning measures described do not replace the respective regulations valid for the establishment!
- Observe the instructions for use of the manufacturers of the disinfectants, particularly the information regarding concentration and material compatibility.
- Surface disinfection of the device should only be performed with the agents specified by the manufacturer, in order to avoid damage to the device.

- **CAUTION!**
Some disinfectant solutions stain the lid components and connecting tubes somewhat. This has no impact on the properties of the materials though.
- Basically, the reprocessing instructions given in Section 7 must be observed when the device is used for another patient.
- The cleaning and maintenance instructions are the result of long years of experience. Shorter cycles may be required, depending on the application and the users' experiences.
- The inhalation device and the accessories must be stored under hygienic conditions!

CAUTION!

- **The device must be disconnected from the mains for cleaning!**
- **The instructions specified in the instructions for use must be observed, particularly as regards the recommended agents.**

6.1.1 AIR FILTER

The air filter (made of foam material) cleans the inhaled air, thus protecting the user and the device.

- The air filter medium cannot be cleaned or disinfected
- For reasons of hygiene, we recommend replacement at least every 4 weeks. If the filter is strongly contaminated, it must be replaced immediately. If the air filter is moist, or if you detect deposits, please observe Section 7.1 of these instructions. Always use the original air filter. Always use an original air filter.

CAUTION!

- **The NEBUJUNIOR® Penguin must not be operated without an air filter.**
- **Make sure you always have a sufficient number of spare air filters available.**

Figure 21

Insertion of the air filter:
Carefully lever out the air filter insert from below (Figure 21).

Figure 22

Place the air filter correctly into the air filter insert (Figure 22). Close the air filter insert by pressing the air filter cover carefully into the filter insert. Take care to insert the air filter properly.

6.0 CLEANING, DISINFECTION, REPROCESSING

EN

6.1.2 AIR TUBE

- If condensation is visible in the air tube, disconnect it from the nebulizer after treatment. The air tube must remain connected to the device.
- Switch on the device. Allow it to operate until the air circulating through the tube has removed any moisture (condensation).
- The air tube must be replaced at least every 6 months (risk of microbial contamination!).
- In the case of contamination, discolouration or damage, replace the air tube immediately (risk of infection!).

6.1.3 NEBULIZER AND ATTACHMENTS

The nebulizer and the attachments used must be cleaned and disinfected immediately after every inhalation routine.

Be careful when disassembling the parts of the nebulizer and avoid damage.

- Remove the mouthpiece, the adapter with the nosepiece, the inhalation mask or the tracheostoma mask, as well as the PEP cap as appropriate, from the nebulizer.
- Disconnect the air tube from the nebulizer and from the air connection of the device.
- Pull off the top section of the nebulizer by twisting both parts (anti-clockwise) gently against each other.
- Remove the blue nebulizer nozzle on the inside.
- Remove any residual medication.
- Under warm running water from the tap, clean the individual parts of the nebulizer as well as the attachments used.
- Check whether dried residue is present. If dried residue is present, place the individual parts in warm water and repeat the previous step.
- To improve the cleaning effect, commercially available dishwashing detergent can be added to the warm rinsing water. If the parts are strongly contaminated, they should be placed in the rinsing water for a longer period of time. Tenacious contaminations must be removed with a soft brush or a soft cloth.
- Place the individual parts in a disinfectant (see the instructions provided with the disinfectant solution for information on the application time).
- Check for any residual disinfectant and remove it.
- Dry the cleaned individual components on a clean base surface. Residual moisture can increase the risk of germs accumulating.
- Alternatively, you can also clean nebulizer incl. PEP cap, mouthpiece, adapter, nosepiece and inhalation mask (exception: tracheostoma mask for tracheotomised patients) in your washing machine (do not clean together with dirty dishes!). Choose a program that runs at a minimum of 55°C for this purpose.
- The parts can also be boiled in hot water at 100°C.
- We recommend cleaning the mouthpiece, the adapter, the nosepiece of the inhalation mask, the nebulizer, as well as the PEP cap after each use.
- Do not use PEP cap, mouthpiece, adapter, nosepiece, inhalation mask and tracheostoma mask if they are damaged or if there are residues/soiling that cannot be removed by rinsing.

6.1.4 DEVICE SURFACE

- In the event of contamination, but at least once a week, the entire surface of the device should be cleaned with a moist cloth (never use a wet cloth) and should then be disinfected with one of the surface disinfectants listed on page 14.
- Some disinfectants can cause discolouration of the material, which does not affect the functionality of the device though.
- Make sure that the device has first been disconnected from the mains power supply!
- The device must never be autoclaved, rinsed under running water or immersed in liquids!
- Ensure that no liquid enters the ventilation slits and reaches the device switch.
- If liquid has penetrated into the device, it may only be started up again after it has been checked by the manufacturer.

6.0 CLEANING, DISINFECTION, REPROCESSING

6.2 DISINFECTION

The different device components can be disinfected with specific disinfectants:

- Nebulizer (instrument disinfectant),
- Attachments (instrument disinfectant),
- Device surface (surface disinfectant).

In this regard, observe the instructions already mentioned in this section as well as the respective disinfectant manufacturers' instructions for use.

Thorough cleaning should always be performed prior to disinfection, if such is required. As a general rule, disinfection is advisable to prevent cross-infection and is also advisable if the device is used in inpatient situations (e.g. hospitals, nursing homes and/or other health care facilities), in order to limit the risk of infections.

6.3 RECOMMENDED INSTRUMENT DISINFECTANTS

Disinfectant	Ingredients	(in 100 g)	Manufacturer
GIGASEPT FF neu (application concentrate)	Succinic acid dialdehyde Dimethoxy tetrahydrofurane Phosphonates, anionic tensides, non-ionic tensides, scents, methyl isothiazolinone	11.9 g 3.2 g	Schülke & Mayr GmbH, Norderstedt
Sekusept Plus (application concentrate)	Glucoprotamin	25.0 g	Ecolab Deutschland GmbH, Düsseldorf

6.4 RECOMMENDED SURFACE DISINFECTANTS

Disinfectant	Ingredients	(in 100 g)	Manufacturer
Dismozon pur steril (granulate)	Magnesium monoperoxyphthalate hexahydrate	80 g	BODE Chemie GmbH, Hamburg
Kohrsolin FF (application solution)	Glutaral Benzyl-C12-18-alkyldimethyl ammonium chlorides Didecyl dimethyl ammonium chloride	5 g 3 g 3 g	BODE Chemie GmbH, Hamburg
Mikrozyd sensitive wipes (ready to use disinfecting wipes)	Alkyl(C12-16)dimethylbenzylammoniumchlorid (ADBAC/BKC (C12-16)); Didecyl-dimethylammoniumchlorid (DDAC), Alkyl(C12-14)ethylbenzylammoniumchlorid (ADEBAC (C12-C14))	0,26 g 0,26 g 0,26 g	Schülke & Mayr GmbH, Norderstedt
Perform (concentrate)	Pentapotassium bis(peroxymonosulphate) bis(sulphate) Anionic tensides, non-ionic tensides, phosphonates, soap, scents	45 g	Schülke & Mayr GmbH, Norderstedt
Bacilloi 30 Foam (application solution)	Propane-1-ol Propane-2-ol Ethanol	45 g 25 g 4.7 g	BODE Chemie GmbH, Hamburg
Mikrobac forte (concentrate)	Benzyl-C12-18-alkyldimethyl ammonium chlorides N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine	19.9 g 5 g	BODE Chemie GmbH, Hamburg

If disinfectants containing aldehyde and disinfectants containing amines are used on the same object, this can cause discolourations.

7.0 MAINTENANCE / REPROCESSING FOR REUSE IN ANOTHER PATIENT

7.1 MAINTENANCE / REPROCESSING FOR REUSE

The NEBUJUNIOR® Penguin is suitable for use in the home as well as in facilities such as clinics, physicians' offices, etc.

To ensure safe and proper operation, a visual and functional check must be performed before each use of the device to determine any possible damage. The cleaning, disinfection and care instructions as described in Section 5 are to be observed. The replacement intervals for accessories and wear components must also be observed. The device may only be used if it is free of defects. There are no additional manufacturer instructions for maintenance.

In the commercial field, checks performed on a regular basis may be required in the scope of accident prevention regulations that apply at the workplace. For non-stationary electrical equipment, these checks are due at least every 24 months. Information in this regard can be obtained from the competent occupational insurance associations.

If the device is to be used for a new patient, it must be professionally reprocessed by the manufacturer or an authorised specialist company in accordance with the relevant medical device standards. Reprocessing for reuse requires the performance of safety-related functional checks by qualified personnel.

For your partners and customers, Andreas Fahl Medizintechnik-Vertrieb GmbH offers unproblematic and fast reprocessing of devices.

- Maintenance / repair of the device, as well as the reprocessing for reuse, may only be performed by the manufacturer or an authorised specialist company.
- This device can be returned for maintenance/repair to Andreas Fahl Medizintechnik-Vertrieb GmbH either directly or via the dealer from whom you have purchased the device.
- Remove all accessory parts before shipping the device. The device itself must be subjected to a surface disinfection.

Liability of Andreas Fahl Medizintechnik-Vertrieb GmbH for personal injury and property damage is excluded if, among other things:

- no original parts are used,
- the instructions regarding proper use in these instructions for use are disregarded,
- assembly, readjustments, modifications, extensions and repairs are performed by specialists who have not been authorised by Andreas Fahl Medizintechnik-Vertrieb GmbH.
- ENo warranty claims are accepted in the event of damages or malfunctions that have arisen due to the use of accessories or consumables from external companies.
- The regulations and instructions that apply to the respective field of application must be observed.

7.2 REPROCESSING FOR REUSE

The handling of the inhalation device has a decisive impact on the reliability and safety of the device. The hygiene measures described in the preceding sections are necessary to protect patients and users and to maintain the functional reliability of the inhalation device. They do not replace reprocessing by the manufacturer.

If the device is to be used for another patient, professional reprocessing at the manufacturer's establishment or at an authorised specialist company in accordance with the Medical Devices Operator Ordinance (MPBetreibV), the Medical Devices Act (MPG) and the Federal Association for Medical Technology (BV-Med) guideline is required.

7.2.1 AIR FILTER

- The air filter must be replaced each time when the device is used for another patient.

7.2.2 NEBULIZER

- The nebulizer, including the PEP cap, must be replaced each time when the device is used for another patient.

7.2.3 AIR TUBE

- The air tube must be replaced each time when the device is used for another patient.

7.2.4 ATTACHMENTS AND ADAPTERS

- The mouthpiece, nosepiece, inhalation mask, tracheostoma mask and adapter must be replaced each time when the device is used for another patient.

7.2.5 DEVICE SURFACE

- If the device is to be used for another patient, the entire surface of the device must be cleaned with a moist cloth (never use a wet cloth) and must then be disinfected with one of the surface disinfectants recommended on page 14.

Some disinfectants can cause discolouration of the material, which does not affect the functionality of the device though.

Make sure that the device has first been disconnected from the mains power supply!

The device must never be autoclaved, rinsed under running water or immersed in liquids!


8.0 TROUBLESHOOTING

EN

The NEBUJUNIOR® Penguin was subjected to a thorough inspection prior to delivery. Should a malfunction nevertheless arise, you may possibly be able to rectify it yourself if you observe the following instructions.

Symptom of malfunction	Possible cause	Remedy
Device fails to start operating	The mains plug is not well seated. Motor has overheated as a result of insufficient cooling.	Check all plug connections. Observe the indicator light; this must be lit if connection has been correctly established. Switch off the device, disconnect the mains plug! Keep the ventilation slits free from obstructions. Wait until the motor has cooled down.
Medication is not sufficiently nebulized	Air tube is kinked. The nebulizer nozzle is blocked. Fluid in air tube. Quantity of medication is too low.	Use the air tube in such a way that no kinks can occur. Clean in accordance with Section 6.1.3/ 6.1.4 or replace the nebulizer. Please see 6.1.2 Check dosage.
Air filter cannot be inserted.	The air filter has not been inserted correctly.	Rotate the air filter through 180 degrees and reinsert the air filter into the opening.

9.0 ACCESSORIES, CONSUMABLES

	Nebulizer REF: 50000-09
	Air tube REF: 50000-01
	Air filter (pack of 10) REF: 50400-07
	Inhalation mask for adults REF: 50000-05
	Inhalation mask for children REF: 50000-08
	Nosepiece REF: 50000-06
	Mouthpiece REF: 50000-10
	Adapter for 15 mm-connector (UNI) REF: 50000-11
	PEP cap
	Carrying bag
OPTIONAL ACCESSORIES (not included in the scope of delivery)	
	LARYVOX® mask Tracheostoma mask for tracheotomised patients REF: 75200

10.0 TECHNICAL DATA

EN

Housing dimensions	300 x 180 x 115mm (11.8" x 7.1" x 4.5")
Weight	2.1kg (4.6Lbs)
Compressor pressure range	30 to 36 psi (2.1 to 2.5 bar)
Compressor working pressure	8 to 16 psi (0.5 to 1.0 bar)
Airflow	6 ~ 8 l/min
Aerosol droplet size	0.5 to 10 µm
Droplet size according to MMAD	Less than 4 µm
Maximum medication capacity	6 ml
Nebulisation capacity	>0.5 ml/min in the case of continuous inhalation (locking device pressed or secured)
Operating noise	Approx. 60 dBA
Ambient conditions	
Transport/storage	-20 C to 70 C (-4 F to 158 F) 10 to 95 % relative humidity
Operation	10 C to 40 C (50 F to 104 F) 30 to 75 % relative humidity
Input voltage	AC 230 V, 50 Hz
Output current	0.8 A
Power admission	180 VA
Classification according to EN 60601-1	Class II (double insulated) Type BF Permanent operation
Protection type	IP 21

11.0 DISPOSAL

EN

- The NEBUJUNIOR® Penguin does not comprise any dangerous goods.
- The housing material is fully recyclable.
- Device and accessories must be decontaminated prior to disposal, since pathogen-containing secretion residues can lead to a hazard.
- Observe careful separation of materials.
- Observe country-specific disposal regulations (e.g. waste incineration).

Disposal in the EU

The inhalation device described above is a high-quality medical device with a long service life. After the end of a life cycle, the device must be transferred for technically correct disposal. According to the EU directives (WEEE and RoHS), the device must not be disposed of via the general household waste. Please observe the existing laws and regulations for disposal of old equipment which apply in the respective country.

Disposal in Germany

In the Federal Republic of Germany, the procedure for disposal of old electric equipment is specified in the Electrical and Electronic Equipment Act (ElektroG). The NEBUJUNIOR® Penguin is exempt from the Electrical and Electronic Equipment Act (ElektroG) in accordance with the EAR (Waste Electrical Equipment Register Foundation) regulation for this type of device. In order to ensure that the regulation-compliant disposal procedure is observed, either return the old device to your responsible specialist dealer, or send it directly for technically correct disposal to Andreas Fahl Medizintechnik-Vertrieb GmbH.

Prior to disposal or transport, all accessory and tube components must be thoroughly cleaned and disinfected. The device itself must be subjected to a surface disinfection.

12.0 INSTRUCTIONS WITH REGARD TO EMC

12.1 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

Der NEBUJUNIOR® Penguin is intended for operation in the electromagnetic environment specified below. The customer or user of the NEBUJUNIOR® Penguin should ensure that the device is used in such an environment.

Electromagnetic immunity tests	IEC 60601- test level	Accordance level	Electromagnetic environment - guidelines
Conducted disturbances as per IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	10 V	<p>Portable and mobile radio devices should not be used at a distance closer to the NEBUJUNIOR® Penguin device, including the cables, than the recommended safety distance as calculated with the equation that applies to the transmission frequency.</p> <p>Recommended safety distance:</p> <p>$d = 0.35 \sqrt{P}$</p> <p>$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the nominal power of the transmitter in watt (W) according to the specifications of the manufacturer of the transmitter and d is the recommended safety distance in metres (m).</p> <p>The field strength of stationary radio transmitters should be lower than the accordance level (b) for all frequencies as determined by an on-site investigation (a).</p> <p>Interference is possible in the vicinity of devices marked with the following warning sign.</p> 
Radiated HF disturbances as per IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	

<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not be applicable in all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.</p>	
a	<p>The field strength of stationary transmitters such as base stations of wireless phones and land mobile radio units, amateur radio stations, AM and FM radio and television transmitters, cannot be predetermined exactly in theory. To establish the electromagnetic environment in regard to stationary transmitters, a study of the location should be considered. If the field strength measured at the location where the NEBUJUNIOR® Penguin is used exceeds the above accordance levels, the NEBUJUNIOR® Penguin should be observed in order to verify that the device functions as intended. If unusual performance characteristics are observed, additional measures may be required, such as for instance a change of orientation or a relocation of the NEBUJUNIOR® Penguin.</p>
b	<p>Above the frequency range of 150 kHz to 80 MHz, the field strength should be below 3 V/m.</p>

12.0 INSTRUCTIONS WITH REGARD TO EMC

EN

The NEBUJUNIOR® Penguin must not be used in an arrangement where it is placed directly next to other devices or stacked with other devices. If the device has to be operated close to or stacked with other devices, the device should be observed in order to verify its proper use when used in this arrangement.

12.2 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY FOR THE NEBUJUNIOR® PENGUIN

Der NEBUJUNIOR® Penguin is intended for operation in the electromagnetic environment specified below. The customer or user of the NEBUJUNIOR® Penguin should ensure that the device is used in such an environment.

Electromagnetic immunity tests	IEC 60601- test level	Accordance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood or concrete or tiled with ceramic tiles. If the floor is covered with synthetic material, the relative air humidity must be at least 30 %.
Fast transient electric disturbances/bursts as per IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input and output lines	± 2 kV for power supply lines	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment.
Surge voltages as per IEC 61000-4-5	± 1 kV common mode voltage ± 2 kV series mode voltage	± 1 kV common mode voltage	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment.
Voltage dips, short interruptions and variations of the power supply voltage as per IEC 61000-4-11	< 5 % U_T (> 95 % dip of U_T for 0.5 cycle) 40 % U_T (60 % dip of U_T for 5 cycles) 70 % U_T (30 % dip of U_T for 25 cycles) < 5 % U_T (> 95 % dip of U_T for 5 s)	< 5 % U_T (> 95 % dip of U_T for 0.5 cycle) 40 % U_T (60% dip of U_T for 5 cycles) 70 % U_T (30 % dip of U_T for 25 cycles) < 5 % U_T (> 95 % dip of U_T for 5 s)	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment. If the NEBUJUNIOR® Penguin user requires continued functions even when the power supply is interrupted, it is recommended to power the NEBUJUNIOR® Penguin from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	30 A/m	Magnetic fields at the mains power supply frequency should be equivalent to the typical values found in the business and hospital environment.

NOTE U_T is the alternating mains power supply voltage before application of the test levels.

12.0 INSTRUCTIONS WITH REGARD TO EMC

EN

- Medical electrical devices are subject to special precautions with regard to EMC and must be installed in accordance with the EMC instructions specified below.
- Portable and mobile HF communication equipment can influence medical electrical devices.
- The use of accessories, transducers or cables other than those specified can lead to increased emission or a reduced electromagnetic immunity of the device or system.

12.3 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

Der NEBUJUNIOR® Penguin is intended for operation in an environment such as the one specified below. The customer or user of the NEBUJUNIOR® Penguin should ensure that the device is used in an environment of this kind.

Interference emission measurements	Accordance level	Electromagnetic environment - guideline
HF emissions as per CISPR 11	Group 1	The NEBUJUNIOR® Penguin uses HF energy exclusively for its internal function. Therefore, its HF emission is very low and it is unlikely that neighbouring electronic equipment will be disrupted.
HF emissions as per CISPR 11	Class B	The NEBUJUNIOR® Penguin is suitable for use in all facilities, including facilities in living quarters and facilities which are directly connected to a public power supply that also supplies buildings used for residential purposes.
Emission of harmonics as per IEC 61000-3-2	Not applicable	
Emission of voltage fluctuations/ flicker as per IEC 61000-3-3	Not applicable	

12.4 RECOMMENDED SAFETY DISTANCES BETWEEN PORTABLE AND MOBILE HF TELECOMMUNICATION DEVICES AND THE NEBUJUNIOR® PENGUIN

Recommended safety distances between portable and mobile HF telecommunication devices and the NEBUJUNIOR® PENGUIN			
The NEBUJUNIOR® Penguin is intended for operation in an electromagnetic environment with controlled HF disturbances. The customer or the user of the NEBUJUNIOR® Penguin can help avoid electromagnetic interference by keeping the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the NEBUJUNIOR® Penguin – depending on the output power of the communication device as specified below.			
Safety distance in metres depending on the transmission frequency			
Nominal power of the transmitter W	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = [0.35] \sqrt{P}$	$d = [0.35] \sqrt{P}$	$d = [0.70] \sqrt{P}$
0.01	0.035	0.035	0.07
0.1	0.11	0.4	0.22
1	0.35	0.35	0.7
10	1.1	1.1	2.2
100	3.5	3.5	7.0
For transmitters of which the maximum nominal power is not indicated in the table above, the recommended safety distance d in metres (m) can be determined using the equation that belongs to the respective column, P being the maximum nominal power of the transmitter in watt (W) as specified by the manufacturer.			
NOTE 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2	These guidelines may not be applicable in all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.		

13.0 LEGAL NOTICES

EN

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the NEBUJUNIOR® Penguin or by repairs, if these modifications or repairs were not carried out by the manufacturer or by an authorised specialist company. This applies both to damages to the NEBUJUNIOR® Penguin itself caused thereby and to any consequential damages caused thereby.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

If the NEBUJUNIOR® Penguin is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH shall, as far as legally permissible, be free of any liability, including liability for defects.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

NEBUJUNIOR® is a registered trademark in Germany of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne.

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NEBUJUNIOR® is een in Duitsland gedeponeerd handelsmerk van Andreas Fahl Medizintechnik-Vertrieb GmbH, Keulen.