
LARYVOX® HME/EXTRA HME/O2 CONNECT

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PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Do not reuse
	Medical devices
	CE marking
	CE marking identification number of the notified body
	Manufacturer

LARYVOX® HME / ARTIFICIAL NOSE

1. FOREWORD

These instructions for use are valid for **LARYVOX® HME / LARYVOX® EXTRA HME**.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

2. SAFETY INFORMATION

en

All HMEs are intended for single patient use only and must not be cleaned or disinfected.

HMEs must under no circumstances be inserted directly into the tracheostoma.

If the HME is impermeable (due to secretion) or shows signs of impermeability, it must be removed and disposed of immediately.

Care must be taken to ensure that the airflow is not blocked by clothing or the like (risk of respiratory distress) to enable unhindered breathing.

If the filter material has slipped in the housing, the product must be replaced and disposed of.

The HMEs must not be fixed to the tracheostomy tube while using a nebuliser or an atomiser.

Do not fill any water into the HME.

When used in combination with tracheostomy tubes and base plates, compatibility and a secure fit must be checked with utmost care.

HMEs must not be disassembled and then reassembled as this will impair their function.

Do not use on patients with too low a respiratory volume!

Respiratory resistance and additional dead space can lead to a build-up of carbon dioxide and cause acute respiratory distress.

Do not grease connections, tubes or other accessories in order to avoid fire hazard.

Smoking during oxygen therapy is dangerous and can lead to severe fire-induced injury.

Only use lotions or ointments labelled as oxygen-compatible to avoid the risk of fire or burns.

There is a fire hazard in connection with oxygen devices and oxygen therapies.

Do not use in the vicinity of sparks or open flames.

Patients who have not previously used HMEs may experience increased respiratory resistance during initial use.

It is recommended to increase the application period slowly.

REF: 49860, 49861, 49862, 49863

The speaking button may only be used exclusively by patients with a partially or fully preserved larynx as well as persons wearing voice prostheses (shunt valve)!

3. INTENDED PURPOSE

HMEs serve to humidify the airways and warm the air one breathes.

They reduce the formation of viscous secretions in the lungs and prevent the trachea from drying out.

HMEs filter the breathing air, thereby reducing intrusion of particles into the patient's airways.

The use of an HME creates the respiratory resistance that is necessary to maintain physiological lung function.

Variants with an O₂ connection also allow for the supply of additional oxygen.

4. INDICATION

HMEs are suitable for spontaneously breathing patients with a tracheostoma.

5. CONTRAINDICATION

The HME should be removed if respiratory resistance is too high.

The use in patients with severe exsiccosis (dehydration), excessive secretion and water accumulation in the lungs and airways, as well as very viscous secretion and haemoptysis is not permitted.

This product is not suitable for patients with low respiratory volumes as the additional dead space may lead to CO₂ retention.

It must not be used outside the recommended tidal volume range, as increased dead space can lead to CO₂ retention if the tidal volume is too low and to insufficient humidification if the tidal volume is too high.

For persons with impaired consciousness, patients with limited motor function and/or persons requiring a high level of care, and who are not in a position to remove the product by themselves, monitoring by trained specialist personnel is required to detect excessive breathing resistance in time and to initiate the appropriate measures to safeguard the airways.

6. COMPLICATIONS AND ADVERSE REACTIONS

If the particle filter is clogged or if signs of clogging appear, the **LARYVOX® HME / LARYVOX® EXTRA HME** filter cassette must immediately be removed and disposed of.

Physical exertion can lead to an increased effort to breathe and shortness of breath due to respiratory resistance.

In this case, it is recommended to withdraw the HME for a short period of time and to train using the HME long-term to achieve better acceptance in the future.

7. PRODUCT DESCRIPTION

LARYVOX® HME (REF 49800) / LARYVOX® EXTRA HME (REF 49860)

- Normal respiratory resistance
- For patients directly after laryngectomy
- For situations with normal physical exertion

LARYVOX® EXTRA HME MEDIUM (REF 49862)

- Average respiratory resistance
- Physical activities

LARYVOX® HME HighFlow (REF 49810) / LARYVOX® EXTRA HME HighFlow (REF 49861)

- Low respiratory resistance
- For new users of HMEs
- For situations with strong physical exertion

LARYVOX® EXTRA HME SPORT (REF 49863)

- Very low respiratory resistance
- For sports activities such as, for example, endurance sports

LARYVOX® O₂ HME (REF 49802)

- Normal respiratory resistance
- With O₂ connection (the oxygen connecting piece fitted at the side of the housing enables a connecting tube to be connected, thereby allowing combination with an oxygen source).
- Ideal for postoperative application

LARYVOX® O₂ CONNECT (REF 49801)

- Normal respiratory resistance
- With O₂ connection (the oxygen connecting piece fitted at the side of the housing enables a connecting tube to be connected, thereby allowing combination with an oxygen source).

LARYVOX® HME / LARYVOX® EXTRA HME filter cassettes feature:

- Plastic housing with an outer diameter of 22 mm
- Compatible with all commonly used 22 mm receptacle systems.
- A foam filter for heat and moisture exchange (see picture 1a)

LARYVOX® EXTRA HME filter cassettes feature:

- A speaking button (see picture 1b), by means of which the tracheostoma can easily be closed for speaking
- A finger groove (see picture 2a) to enable easy removal of the HME filter cassette from the tracheostoma housing for coughing
- Side openings for the airflow to enable unhindered and easy breathing (see picture 2b).

8. APPLICATION

en Check whether the packaging and the product are in a faultless condition.

Check the use-by or expiry date. Do not use after this date.

Do not use the product after this date.

Check the stated tidal volume.

Fit (see picture 4) the **LARYVOX® HME / LARYVOX® EXTRA HME / LARYVOX® O₂ CONNECT** filter cassettes vertically onto the 22-mm connection of the tracheostomy tube or base plate, exerting gentle pressure.

Check the respiratory resistance.

Remove the **LARYVOX® HME / LARYVOX® EXTRA HME / LARYVOX® O₂ CONNECT** filter cassettes after a (see picture 7) wearing period of 24 hours at most, or if the filter cassette is contaminated with secretions for instance.

To do this, hold the base plate or the tracheostomy tube in position with two fingers (see picture 8) and remove the filter cassette by tilting it downward slightly.

LARYVOX® EXTRA HME

Press the speaking button down in order to check whether it is fully functional (see picture 3a). When the speaking button is released, it immediately returns to the open starting position (see picture 3b).

For speaking, press the speaking button (see picture 5) down slightly with your finger so as to enable phonation by means of the closure.

When the speaking button is released, it immediately returns to its open breathing position (starting position) (see picture 6).

9. HYGIENE INSTRUCTIONS

HMEs are single-patient products and intended for single use only.

They must not be cleaned or disinfected, as use of disinfectant can render the HME unusable.

In addition, bacterial or fungicidal colonisation of the filter material can lead to an increased risk of infection.

Cleaning, disinfection or (re-)sterilisation as well as use longer than 24 hours can impair the safety and function of the product and are therefore not permissible!

10. STORAGE

Store in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

The maximum period of use / wearing period is 24 hours.

The filter must be changed even several times a day if required, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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