




LARYVOX® PAD

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I. FOREWORD

These instructions for use are valid for LARYVOX® PAD. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

II. INTENDED USE

LARYVOX® PAD is used to align the skin level of an irregular tracheostoma that is situated below the skin level.

The speech option is improved by reducing undesired air loss.

The product is intended for domestic use only.

III. WARNINGS

LARYVOX® PAD is a single-patient product and is thus intended exclusively for one-time or repeated use in one individual patient.

When skin irritations occur, the LARYVOX® PAD must be changed immediately.

When LARYVOX® PAD is used for the first time, its suitability for sealing the tracheostoma should be checked.

In problematic individual cases, sealing of the tracheostoma may not be sufficient.

LARYVOX® PAD is intended exclusively for use in combination with LARYVOX® Tapes.

If you wish to use LARYVOX® PAD in combination with a tracheostomy tube, please discuss this with your attending physician beforehand.

IV. CONTRAINDICATIONS

Do not use LARYVOX® PAD in case of allergy or intolerance to the material.

Do not wear LARYVOX® PAD during laser treatment or electrocautery.

There is a considerable risk of injury from exposure to heat and, if applicable, from the release of toxic combustion products.

Use of the LARYVOX® PAD should be avoided while undergoing radiotherapy due to the risk of skin irritations.

After undergoing radiotherapy, only use LARYVOX® PAD after consulting your treating physician. Start using the tape at the earliest 2 weeks after completing radiotherapy.

V. PRODUCT DESCRIPTION

LARYVOX® PAD is a sealing plate made of soft, skin-coloured medical silicone.

It is used to align the skin level of an irregular tracheostoma that is situated below the skin level.

LARYVOX® PAD can also be used as a seal between the tracheostoma and tracheostomy tube.

Please consult your physician on this matter!

Due to its material properties and ergonomic shape, the LARYVOX® PAD adapts to the skin around the tracheostoma as well as to the shield of the tracheostomy tube used.

LARYVOX® PAD is particularly skin-friendly due to the material used.

No adhesive is required to fix the pad.

The plano-convex shape of the LARYVOX® PAD ensures high quality sealing and at the same time provides cushioning between the tracheostomy tube and the tracheostoma.

In addition to reduced leakage of secretion, the lateral flow of breathing air can also be minimised.

At the same time, the tracheostomy tube can be stabilised in the tracheostoma.

This makes LARYVOX® PAD suitable for ventilated patients as well as for improving phonation, e.g. when wearing a shunt valve.

VI. INSTRUCTIONS

1. The skin surrounding the tracheostoma must be clean, dry and free from grease.
2. Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or OPTIFAHL® stoma cleansing wipes (REF 33200).
3. Dry the skin.
4. When opening the packaging, make sure that the LARYVOX® PAD is undamaged and that the indicated use-by date has not been exceeded.
5. After removing the protective foil backing ❶, place the LARYVOX® PAD with the adhesive side centred over the tracheostoma ❷.
6. If you use LARYVOX® Pad in combination with a tracheostomy tube, you can also push the cannula shaft through the opening of the pad before insertion (❶a).
7. Otherwise, insert the cannula as usual after fixing the LARYVOX® PAD.
8. Finally, place LARYVOX® Tape over the LARYVOX® PAD according to the instructions for use ❸.
9. For the best possible hold, use a tape with a sufficiently large adhesive area (REF 48200, 48300, 48220 48320, 48230, 48330, 48240, 48340, 48150, 48250, 48400).
10. When fixing an inserted tracheostomy tube or a button with a cannula strap, please observe the respective instructions for use.
11. Remove the LARYVOX® PAD and LARYVOX® Tape by carefully removing the adhesive strips from the skin ❹.

CAUTION!

Do not use any fatty creams!

Otherwise, secure attachment of the Pad and Tape is no longer ensured.

VII. CLEANING

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. Please consult your physician.

Cleaning Steps

Rinse the LARYVOX® PAD thoroughly under lukewarm running water.

Then allow LARYVOX® PAD to dry in a clean environment.

If you are using LARYVOX® PADs with a tracheostomy tube, the Pad should be removed and cleaned each time the tube is changed.

VIII. STORAGE

Currently not used and cleaned LARYVOX® PADs should be stored in a clean plastic container in a dry environment away from dust, sunlight and/or heat.

Still packaged replacement pads should be stored in a dry environment away from sunlight and/or heat.

We highly recommend keeping a replacement pad on hand to ensure continuous supply.

IX. PERIOD OF USE

After 29 days, the LARYVOX® PAD should be disposed of properly and replaced with a new product.

Deviations from the expected wearing period occur as a result of individual wearing conditions and skin properties.

Should the material show cracks or other damage, the product must be replaced immediately.

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Product specifications described herein are subject to change without notice. LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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