# NebuJunior® Bär Inhalator

REF 50300

CE 0482

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## Operating instructions

Please read the instruction manual carefully before initial operation. Keep instruction manual for further use!

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<b>12</b>	<b><i>Safety Precautions</i></b> .....	<b>45</b>
12.1	<b>Intended use</b> .....	<b>45</b>
12.2	<b>Electrical connection / Failure</b> .....	<b>47</b>
12.3	<b>Device inspection</b> .....	<b>47</b>
12.4	<b>Device modifications and spare parts</b> .....	<b>47</b>
12.5	<b>Product enhancements</b> .....	<b>48</b>
<b>13</b>	<b><i>Special Safety Precautions</i></b> .....	<b>49</b>
13.1	<b>Device-specific hazards and residual risks</b> .....	<b>49</b>
<b>14</b>	<b><i>Device Description</i></b> .....	<b>50</b>
14.1	<b>Device overview / Scope of supply</b> .....	<b>51</b>
14.2	<b>Optional accessories (not included)</b> .....	<b>53</b>
14.3	<b>Inhalation device Components and their Functions</b> .....	<b>54</b>
14.4	<b>Product Specifications</b> .....	<b>55</b>
<b>15</b>	<b><i>Operation of the NebuJunior<sup>®</sup> Inhalation device</i></b> .....	<b>56</b>
15.1	<b>General Requirements</b> .....	<b>56</b>
15.2	<b>Initial operation</b> .....	<b>56</b>
15.3	<b>Inhalation Therapy</b> .....	<b>57</b>
15.3.1	Preparation of the inhalation therapy .....	57
15.3.2	Starting the inhalation therapy.....	61
<b>16</b>	<b><i>Cleaning and Disinfection</i></b> .....	<b>63</b>
16.1	<b>Preparations for the cleaning process</b> .....	<b>64</b>
16.2	<b>Cleaning and disinfection after the inhalation therapy</b> .....	<b>64</b>
16.3	<b>Care and maintenance of the air tube</b> .....	<b>68</b>
16.4	<b>Exchanging the air filter</b> .....	<b>69</b>
16.5	<b>Recommended disinfectants</b> .....	<b>70</b>
16.6	<b>Device inspection after cleaning/disinfection</b> .....	<b>71</b>

16.7	Storage .....	71
17	<i>Technical Faults/Malfunctions and Troubleshooting</i> .....	72
17.1	Malfunctions to be corrected by the user/patient.....	72
17.2	Technical faults only to be corrected by the Manufacturer .....	73
18	<i>Spare Parts / Accessories for NebuJunior®</i> .....	74
19	<i>Maintenance / Processing for Re-use</i> .....	75
20	<i>Disposal</i> .....	76
21	<i>Declaration of conformity</i> .....	78
22	<i>Guarantee Certificate</i> .....	79

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CE 0482

## Operating instructions

**Please read the instruction manual carefully before initial operation. Keep instruction manual for further use!**

The NebuJunior® inhalation device is a modern inhalation device used for respiratory air humidification, which operates by compressed air. All types of medication permissible for an inhalation therapy can be nebulised with the NebuJunior® inhalation device. Please consult your physician for the suitable medication.

Inhalation is possible either with a mouthpiece or with an inhalation mask - both accessories are included. A special adapter tube with tracheostoma mask is available for tracheostomised patients.

For cleaning and maintenance, the **NebuJunior®** inhalation device can easily be disassembled.

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## **12 Safety Precautions**

### **12.1 Intended use**

The operational safety of the NebuJunior® inhalation device can only be guaranteed if the device is used appropriately (intended use).

**The NebuJunior® inhalation device is only used appropriately if:**

- the NebuJunior® is used in accordance with these instructions, for the prophylaxis and treatment of illnesses of the upper and lower respiratory tract, e.g. for domestic inhalation therapy or in medical facilities; AND
- there is no medical contraindication for an inhalation therapy – if necessary, please consult your physician for advice; AND
- the medication used and its dosage were determined by your physician; AND
- children and patients in need of medical assistance use the inhalation device only under supervision.

**Please carefully observe the operating and maintenance instructions!**

**Any inappropriate use of the NebuJunior® inhalation device is not permitted and includes:**

- operating the NebuJunior® in a damaged condition;
- inhaling medication and remedies that are not permissible for inhalation therapies;
- technical modifications of the NebuJunior® inhalation device;
- the use of the NebuJunior® inhalation device by unattended children;
- the use of the NebuJunior® Inhalator in damp rooms or the shower / bathtub;
- the outdoor use of the NebuJunior® inhalation device;
- the Non-observance to any part of the operating instructions;
- the reprocessing of the device for reuse and resale without authorisation of the manufacturer.

**Any inappropriate use can lead to severe material damages and personal injuries!**

## 12.2 Electrical connection / Failure

**In case of failure, the NebuJunior® inhalation device must be completely disconnected from the electrical power by pulling the plug out of the power socket. Please observe the details in Chapter 6 (Technical Faults and Disturbances)!**

## 12.3 Device inspection

Before operating the NebuJunior® inhalation device and after each cleaning and/or disinfection procedure, please make sure that:

- all parts of the device are assembled properly
- that the air tube (4) is not bent or otherwise defect

See also the chapters concerning cleaning, maintenance and special safety advice.

## 12.4 Device modifications and spare parts

**Caution!**  
**Inadequate or inapplicable accessory parts may lead to injuries! Please only use the manufacturer's original accessories and spare parts!**

Any modifications or rebuilding procedures on the NebuJunior® inhalation device on the user's own authority lead to an immediate expiration of the manufacturer's liability and warranty obligations. The same applies for the utilisation of non-original spare parts.

- Repair and maintenance works may only be carried out by the manufacturer himself or by his authorised representatives.

## **12.5 Product enhancements**

We are constantly working on our devices in order to offer our customers only products that correspond to the state-of-the-art technology and to all applicable safety requirements. Therefore, please inform us about any problems with the NebuJunior<sup>®</sup> inhalation device.

## 13 Special Safety Precautions

### 13.1 Device-specific hazards and residual risks

So-called residual risks are dangers which cannot be eliminated during the construction of the device.

**Caution!**

**Risk of infection due to viruses, bacteria, germs and fungi; may occur in case of insufficient cleaning and/or disinfection.**

## 14 Device Description



The aim of this chapter is to make the patient familiar with the display panel, the control elements and the overall construction of the NebuJunior<sup>®</sup> inhalation device.

**Prior to the treatment of any illnesses, we highly recommend a medical examination by a physician. The medication used and its dosage must be determined by your physician.**

1. The NebuJunior<sup>®</sup> inhalation device is used for the prophylaxis and treatment of illnesses of the upper and lower respiratory tract.
2. The medicinal vapour (aerosol) is generated by compressed air.
3. In principle, every kind of medication available and authorised for inhalation therapy may be used with the NebuJunior<sup>®</sup> inhalation device.

### 14.1 Device overview / Scope of supply


**6 Nebuliser**


**2 Air tube**

**8 Inhalation mask (for adults)**

**16 Inhalation mask (for children)**

**13 Air filter**

**5 Mouthpiece**

**9 Nosepiece**

## 14.2 Optional accessories (not included)

**14 Laryvox® mask,  
Tracheostoma mask for  
tracheostomy patients**

**15 Adapter for tracheostoma  
mask**

### 14.3 Inhalation device Components and their Functions

Pos.	Component	Function
1	Air compressor	generates compressed air
2	Air tube	transmits compressed air to nebuliser
3	Air filter compartment	opening on top of device for insertion of the air filter
4	Air tube connection	for attaching the air tube to compressor
5	Mouthpiece	for mouth inhalation
6	Nebuliser	nebulises the inhalant / medication
7	Power cable	for connection to the power socket (230V )
8	Inhalation mask (for adults)	for inhalation treatment via the patient's face
9	Nosepiece	for inhalation treatment via the patient's nose
10	Air vents	vents for air supply
11	Power switch	for switching inhalation device on / off
12	Handle	for carrying compressor
13	Air filter	for filtering/processing the air
14	Tracheostoma mask <b>Optional accessory</b>	for inhalation via the tracheostoma
15	Adapter for Tracheostoma mask <b>Optional accessory</b>	for connecting the tracheostoma mask to the nebuliser
16	Inhalation mask (for children)	for inhalation treatment via the patient's face

## 14.4 Product Specifications

Description	Data
Power supply	AC 230V, 50Hz
Power consumption	80VA +/- 20VA
Max. filling quantity (medication)	13ml
Particle size aerosol	0.5 - 10 µm
Particle size acc. to MMAD	smaller than 4 µm
Sound level	Ca. 60 dBA
Nebulization capacity	0.2ml/min
Pressure range of compressor	29 - 43 Psi (2.0 - 3.0 bar)
Operating Pressure of compressor	8 - 16 Psi (0.7 - 1.1 bar)
Air flow	6 ~ 8 l/min
Operating temperature	10°C - 40°C (50°F - 104°F)
Rel. Air humidity during operation	10 - 95% RH
Storage temperature	-20°C - 70°C (-4°F - 158°F)
Rel. Air humidity during storage	10 - 95 % RH
Dimensions	175 x 120 x 180mm (6.9" x 4.7" x 7.1")
Weight	1.6kg (3.5Lbs)
Standard accessories	Power cable, nebuliser, mouthpiece, air tube, air filter (5 pcs), inhalation mask (for adults), inhalation mask (for children), carrying bag
Optional accessories	Tracheostoma mask, adapter for tracheostoma mask
Device class	Type BF
IP class	II
Operation period	Short term operation – max. 30 minutes
CE mark	CE 0482

Nebuliser, air tube, mouthpiece, nosepiece, inhalation mask (for adults), and inhalation mask (for children), are made of PVC.

Specifications are subject to change without prior notice.

## 15 Operation of the NebuJunior® Inhalation device

This chapter informs you about safety precautions which need to be observed prior to each inhalation therapy.

### 15.1 General Requirements

The NebuJunior® inhalation device is ready for use after the packaging has been removed and after the following instructions have been carried out:

- Inspect the device for any external visible damages.  
The device must be free of damages.
- Do not operate the NebuJunior® inhalation device in case you detect any damages to the device.
- Place the NebuJunior® inhalation device on a flat and stable surface, e.g. a table.
- Connect the NebuJunior® inhalation device to a 230V / 50Hz power socket.
- Do not place the NebuJunior® inhalation device on a tablecloth or any piece of fabric.
- Check that all air vents are unobstructed.

### 15.2 Initial operation

Prior to initial operation, make sure that all components as listed in Chapter 14.1. (Device overview / Scope of delivery ) are included in the NebuJunior® inhalation device.

Then clean and disinfect the device as described in Chapter 16.


**In order to guarantee a faultless operation, the NebuJunior® inhalation device is at no time to be exposed to direct sunlight or to high temperatures!**

## 15.3 Inhalation Therapy

### 15.3.1 Preparation of the inhalation therapy

**Please follow the instructions in correct order.**

- Make sure that the NebuJunior<sup>®</sup> inhalation device and the nebuliser parts are cleaned and disinfected.
- Make sure that the air tube (2) is not bent or otherwise defect.
- Pull off the top part of the nebuliser (6), by twisting both parts gently in opposite directions (anticlockwise). **Image 1**

**Image 1**

...continued

- Fill the correct amount of medication / inhalant, as prescribed by your physician or as described in the medication package insert, and fill it into the lower part of the nebuliser **(6)**. *Image 2a*.
  - Please note that for safety reasons, the maximum filling quantity for the nebuliser is **13 ml**.
- 
- Re-attach the upper part of the nebuliser **(6)** to the lower part.
  - Re-assemble the nebuliser, by carefully twisting both parts (clockwise). Please make sure that both parts are fixed securely. *Image 2b*

*Image 2a*

*Image 2a*

**...continued**

- Connect the air tube (2) to the opening on the bottom of the nebuliser (6).
- Then attach the air tube to the air tube connection (4) of the compressor. **Image 3**

**Image 3**

**...continued**

- now attach one of the following application-specific accessories to the top of the nebuliser **(6)**:  
the mouthpiece (**Image 4a**) or the nosepiece (**Image 4b**) or  
the inhalation mask (**Image 4c**) or the adapter with  
tracheostoma mask (**Image 4d**).

**Image 4a**

**Image 4b**

**Image 4c**

**Image 4d**

### 15.3.2 Starting the inhalation therapy

- Keep the nebuliser **(6)** in a vertical position. The nebuliser is operable up to an angle of 45 degree.
- Switch on the NebuJunior<sup>®</sup> inhalation device by pressing the power switch **(11)**.
- Breathe the inhalation vapour carefully in a few times so that the respiratory tract can accustom itself to the medicinal vapour.
- During the inhalation process, remain in a relaxed sitting position and breathe calmly and evenly.
- If you use the mouthpiece **(5)** for inhalation, please enclose the mouthpiece firmly with your lips.

**Children and/or people in need of medical assistance may only use the inhalation device with an inhalation mask or tracheostoma mask (the latter in case of tracheostomy) and must be supervised by an adult and/or nursing staff.**

**If the tracheostoma mask (optional accessory) is used for inhalation, the corresponding adapter must also be used. Under no circumstances attach the tracheostoma mask directly to the nebuliser!**

- After having breathed in the medicinal vapour (aerosol), hold your breath for a few seconds and then breathe out slowly. Patients inhaling orally should breathe out via their nose and keep their mouth closed.
- Press the inhalation mask or tracheostoma mask (**Image 4c, 4d**) lightly onto your face or tracheostoma and inhale calmly.
- When exhaling, please hold the inhalation mask / tracheostoma mask (**Image 4c, 4d**) at a small distance (fingerbreadth) away from your face or tracheostoma.

**...continued**

- If not otherwise prescribed by your physician, the inhalation therapy should be performed 2-3 times a day for approx. 10-15 minutes.
- If any mucus is loosened by the inhalation therapy, this should be coughed up and spat out. Tracheotomised patients should remove the mucus by means of a suction apparatus and sterile suction catheters.
- Directly after the inhalation therapy, draughty and cold surroundings should be avoided.
- After having finished the inhalation therapy, please switch off the NebuJunior<sup>®</sup> inhalation device at the power switch **(11)**.
- Then clean/disinfect the inhalation device as described in Chapter **16**.

## 16 Cleaning and Disinfection

This chapter instructs/-explains the user/patient how to keep the NebuJunior® inhalation device and its accessories in a fully functional and good hygienic condition.

**Risk of infection due to viruses, bacteria and germs, infections can transmit serious illnesses. After each inhalation therapy, the device must be cleaned and disinfected in order to minimise the risk of infection.**

All parts which were in direct contact with either the medication or the patient's body must be cleaned and disinfected thoroughly immediately after the inhalation therapy. For the cleaning process, the nebuliser needs to be disassembled into its single components.

The air compressor (1) should be cleaned/disinfected on the exterior by using of a soft, clean cloth. After cleaning/disinfection, all components must be dried thoroughly. Any humidity remaining on the components increases the risk of germ development.

The cleaning process is described in detail on the following pages.

## 16.1 Preparations for the cleaning process

**Danger to life by electrical shock!**

**Always disconnect the power cable from the power socket before starting the cleaning/disinfection process.**

Please proceed with caution when disassembling the nebuliser to avoid damaging.

- Remove the mouthpiece **(5)**, the nosepiece **(9)**, the inhalation mask **(8/16)** or the tracheostoma mask **(14)** with adapter **(15)** from the nebuliser **(6)**.
- Remove the air tube **(2)** from the nebuliser and from the air tube connection **(4)** on the compressor **(1)**.
- Pull off the top part of the nebuliser, by twisting both parts gently in opposite directions (anticlockwise).
- Remove all medication/inhalant residues.

## 16.2 Cleaning and disinfection after the inhalation therapy

Please only use disinfectants and cleaning procedures that have been approved of and authorised by the Robert Koch Institute or by the German Society for Hygiene and Microbiology (DGHM).

You will find a list of recommended disinfectants in Chapter **16.5**. For cleaning the single components of the NebuJunior<sup>®</sup> inhalation device, please proceed according to the tables on the following pages:

Cleaning and disinfection for the parts shown here	Parts	Remarks
<ol style="list-style-type: none"> <li>1. Rinse with warm tap water*.</li> <li>2. Check for encrustations and if necessary, soak these in hot water, then repeat Step 1.</li> <li>3. Immerse components in the disinfection solution (for duration, check instructions for use of the solution).</li> <li>4. Rinse component again to remove disinfectant residues.</li> <li>5. Visual inspection concerning disinfectant residues.</li> </ol> <p>Leave to dry on clean cloth (e.g. a flattened towel).</p> <p>Boiling the components with 100° Celsius hot water is also possible.</p>		<p>In case of reprocessing the inhalation device for re-use (different patients!), this component must be exchanged.          Risk of infection!</p>
		<p>In case of reprocessing the inhalation device for re-use (different patients!), this component must be exchanged.          Risk of infection!</p>
		<p>In case of reprocessing the inhalation device for re-use (different patients!), this component must be exchanged.          Risk of infection!</p>
		<p>In case of reprocessing the inhalation device for re-use (different patients!), this component must be exchanged.          Risk of infection!</p>

		In case of reprocessing the inhalation device for re-use (different patients!), these components must be exchanged. Risk of infection!
		
<b>Store all components in a clean, dry and dust-free environment and protect them from direct sunlight!</b>		

**Electrical components may be damaged by humidity!  
 Do not immerse the air compressor (1) in liquid and  
 do not hold the compressor under running water!**

Cleaning and disinfecting the air compressor (1)	Components
<ol style="list-style-type: none"> <li>1. Switch off the device and disconnect the power cable.</li> <li>2. Clean with a soft, damp cloth.</li> <li>3. Wipe-Disinfection of the compressor surface (see list of recommended disinfectants).</li> </ol> <p><b>Caution!</b>  <b>Do not hold under running water or immerse device in any liquid!</b>  <b>Do not spray disinfectant into air vents or onto the power switch!</b></p> <ol style="list-style-type: none"> <li>4. Dry device.</li> </ol>	
<p><b>Store device in a dry and dust-free environment!</b></p>	

### 16.3 Care and maintenance of the air tube

If you detect condensate accumulations within the air tube, then please detach the air tube from the nebuliser after having finished the inhalation therapy. The air tube must however remain connected to the air compressor. Then switch on the air compressor. Leave the compressor turned on until the air flow has cleared all condensate accumulations. Exchange the air tube every 6 months at least in order to prevent germ development!

In case of reprocessing the inhalation device for re-use (different patients!), pollutions, discolourations or damage the air tube must be exchanged. Risk of infection!

**Store the air tube in a dry and dust-free environment and protect it from direct sunlight!**

## 16.4 Exchanging the air filter

The air filter (*Image 5b*), which is located inside the air filter compartment (*Image 5a*) should be exchanged at least once a month (given a normal operating frequency). In case of reprocessing the inhalation device for re-use (different patients!), the air filter must be exchanged immediately (Risk of infection!).

Please examine the air filter in regular intervals (approx. after every 10th operation). If the air filter is clogged or dirty (gray or brown colouring), then it must be exchanged immediately. The same applies in case the air filter gets damp or wet.


**Never operate the device without a filter. The filter must under no circumstances be cleaned and reused! Only use original NebuJunior® air filters. Otherwise, your air compressor may not work correctly and a satisfactory inhalation therapy can not be guaranteed.**

For exchanging the air filter, please proceed as follows: Remove the cover from the air filter compartment (*Image 5a*). Then take the used air filter out and replace it with a new one. Close the air filter compartment, by carefully pressing the cover back onto the compartment.

*Image 5a*

*Image 5b*

## 16.5 Recommended disinfectants

**Please follow the instructions for use and the safety instructions of each recommended disinfectant. The prescribed dosages and disinfection periods must be observed.**

The below mentioned disinfectants are suited for the disinfection procedure and permit a frequent disinfection without damaging the inhalation device or its components.

<b>Disinfectants for exterior disinfection (air compressor)</b>
<b>D 10</b> Manufacturer <b>orochemie</b> Kornwestheim (Germany). <b>PROPANO AF</b> Manufacturer <b>esteer PHARMA GMBH</b> Reilingen (Germany).

<b>Disinfectants for disinfection bath (nebuliser components)</b>
<b>D 10</b> Manufacturer <b>orochemie</b> Kornwestheim (Germany). <b>DESRESON AF</b> Manufacturer <b>esteer PHARMA GMBH</b> Reilingen (Germany).

## **16.6 Device inspection after cleaning/disinfection**

Before each inhalation therapy, please check the functionality of your NebuJunior<sup>®</sup> inhalation device as follows:

- Check all components of the NebuJunior<sup>®</sup> inhalation device for damages or discolourations. Any affected component may have to be exchanged. A list of all spare parts can be found in Chapter **18**.
- Then connect the air tube **(2)** to the air tube connection **(4)** and to the nebuliser **(6)**.
- Turn on the power switch **(11)** on the NebuJunior<sup>®</sup> inhalation device.
- Make sure that air flows out of the nebuliser.

## **16.7 Storage**

For storing the NebuJunior<sup>®</sup> inhalation device, please put the device into the provided carrying bag.

Make sure to store the inhalation device and all its components in a dry and dust-free environment.

## 17 Technical Faults/Malfunctions and Troubleshooting

### 17.1 Malfunctions to be corrected by the user/patient

Disturbance / Fault	Possible Cause	Corrective Action
The NebuJunior® inhalation device cannot be turned on or the engine does not start.	Power Cable not connected.	Make sure the power cord is connected to the power socket.
	Overheating of motor due to too long operating period or due to insufficient cooling.	Switch off inhalation device and disconnect power cable from power socket. Make sure, that the air vents are unobstructed. Wait for motor to cool down.
The medication/inhalant is not sufficiently nebulised.	Air tube is squeezed or bent.	Arrange the air tube so that it is neither squeezed nor bent.
	Nozzle of the nebuliser may be clogged.	Clean nozzle according to chapter <b>16.2</b> or exchange the nebuliser
	Quantity of medication may be too low.	Check dosage.

If the problem still remains after the corrective action was taken, please return the device to Andreas Fahl Medizintechnik-Vertrieb GmbH for inspection and repair works.

**17.2 Technical faults only to be corrected by the Manufacturer**

<b>Disturbance / Fault</b>	<b>Possible Cause</b>	<b>Corrective Action</b>
The NebuJunior® inhalation device cannot be turned on.	Electrical fault.	Return the faulty NebuJunior® inhalation device for repair works to Andreas Fahl Medizintechnik-Vertrieb GmbH.
Device does not generate air pressure.	Sealings / gaskets may be damaged.	For quick and reliable repair of your NebuJunior® inhalation device, please send us a detailed fault description.
Any other disturbance / fault.	Various possible causes.	Please ensure to send the device in a sturdy and secure packaging.

## 18 Spare Parts / Accessories for NebuJunior®

Pos. *	Part description	REF
2	Air tube for Nebufirst®, Nebuport®, NebuJunior®	50000-01
5	Mouthpiece for Nebufirst®, Nebuport®, NebuJunior®	50000-02
6	Nebuliser for Nebufirst®, Nebuport®, NebuJunior®	50000-03
8	Inhalation mask (for adults)	50000-05
9	Nosepiece for Nebufirst®, Nebuport®, NebuJunior®	50000-06
13	Air filter for Nebufirst®, NebuJunior® Bär	50000-07
14	Laryvox® mask, Tracheostoma mask	75200
15	Adapter tracheostoma mask for Nebufirst®, Nebuport®, NebuJunior®	50110
16	Inhalation mask (for children)	50000-08

\* See spare part number according to device overview in Chapter **14.1.**

## 19 Maintenance / Processing for Re-use

The NebuJunior® inhalation device is suited for domestic use in private households and also for health institutions such as clinics or doctor's practices.

If the NebuJunior® inhalation device is used in a professional environment, all applicable regulations for accident prevention must be followed in order to assure the employees' safety. These regulations include an inspection concerning electrical safety, to be effected within the scope of maintenance works by specially trained personnel. The frequency of these maintenance works is not given by the manufacturer. The operator himself must decide on the inspection intervals by taking the duration, the frequency and the mode of operation of the NebuJunior® inhalation device into consideration.

In general, a hygiene-based inspection of the device must be undertaken after every inhalation therapy. If the nebuliser is used by more than one patient it must be cleaned and disinfected before each change of patient. It is absolutely necessary to follow the aforementioned cleaning and disinfection guidelines in order to avoid any risk of infection!

Devices that are being operated solely in a domestic environment and by one patient must only be resold or otherwise transferred to another person on condition of a previous technical revision (complete overhaul), a hygiene-based inspection and maintenance works by the manufacturer **Andreas Fahl Medizintechnik-Vertrieb GmbH** or by one of his authorised representatives. Please ensure that only procedures which are sufficiently validated for the devices and products are used. The same regulation applies for all general repair and maintenance works on the device.

## 20 Disposal

The NebuJunior<sup>®</sup> inhalation device is intended for domestic use in private households (individual user / patient) and in a professional environment (doctor's practice, clinics and hospitals, other health-care facilities).

The device is subject to Category 8 of the German ElektroG (*Act on Waste Electrical and Electronic Equipment*), Paragraph 2, Clause 1, which defines the manufacturer's obligation for recycling and disposal of the concerned electric or electronic equipment.

However, according to Paragraph 2, Clause 1 of the same Act, products that are considered "infectious" are exempt from this regulation.

According to the German ElektroG, devices are considered infectious, if blood or other biologic contaminants (like secretion) are able to infiltrate the device. This regulation affects the tracheostoma mask (**14**) and the adapter (**15**), for which the aforementioned contamination risk cannot be eliminated. Therefore, these accessory parts must be disposed of together with the ordinary household waste.

Please inform the manufacturer if your device is no longer needed. We will advise you on the further procedure concerning disposal / recycling of the NebuJunior<sup>®</sup> inhalation device.

**Prior to disposal of the device, please observe the following:**

**Risk of infection by viruses, bacteria and germs.**

**Infections can transmit critical illnesses.**

**The device must be clean and disinfected prior to disposal in order to minimise the risk of infection to other people.**

### Attention!

Please observe the relevant, national regulations concerning disposal of the device.

In Germany, the regulations according to ElektroG, MPG, etc. apply. Further information can be found on the website **www.stiftung-ear.de**, which is the officially appointed body of the Federal Environmental Agency.

- The NebuJunior<sup>®</sup> inhalation device must be disposed of according to the regulations for electrical waste and **may not** be disposed of with ordinary household waste.
- In Germany, separate waste collection points for electrical devices have been set up by the Public Waste Management Authorities.
- **For other countries, please refer to the national directives and regulations!**

## 21 Declaration of conformity

The Fahl® Inhalation devices are medical devices as defined by the MPG (German medical devices act) and conform to, among others, the "Basic Requirements" as per Annex I of Directive 93/42 EEC.

NebuJunior® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorized product alterations or improper use, care, and/or handling.

Sales and deliveries of all Andreas Fahl Medizintechnik-Vertrieb GmbH products are carried out exclusively according to our General Terms and Conditions of Business which can be received directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

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## 22 Guarantee Certificate

# Guarantee Certificate

By purchasing the NebuJunior® inhalation device, the customer acquired a medical device of first quality. If defects/faults appear within the guarantee period (24 months from purchase date), the manufacturer commits himself to repair or exchange the device free of charge.

Any damages or defects on the device that are caused by improper use or by technical modifications effected by others than the manufacturer or his authorised representatives lead to an immediate expiration of the guarantee.

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Serial no. of device \_\_\_\_\_

Delivery date \_\_\_\_\_

In case of a warranty claim, please return the device together with the completely filled in guarantee certificate to:

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**August-Horch-Straße 4a**  
**51149 Köln/ Cologne (Germany)**

Technical modifications and repair works may only be performed by the manufacturer or his authorised representatives.