











LARYVOX® INSERTER UNIEK PRO

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PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

en

	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Do not reuse
	Medical devices
	Non-sterile
	Do not use if package is damaged
	CE marking
	Manufacturer

LARYVOX® INSERTER UNIEK PRO

1. FOREWORD

These instructions apply to **LARYVOX® INSERTER UNIEK PRO**.

These instructions for use serve to inform the physician and medical professional staff in order to ensure technically correct handling of **LARYVOX® INSERTER UNIEK PRO**.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

The **LARYVOX® INSERTER UNIEK PRO** is a single-patient product and is intended for one single use only.

Never use the **LARYVOX® INSERTER UNIEK PRO** if the device exhibits visible deficiencies or defects.

A damaged voice prosthesis must not be used.

The voice prosthesis must always be placed with care in the shunt channel.

Pay attention to the correct position of the voice prosthesis.

Please also observe the notes / warnings in the instructions for use provided by the voice prosthesis manufacturer.

Do not make more than three placement attempts with a **LARYVOX® INSERTER UNIEK PRO** applicator.

3. INTENDED PURPOSE

INSERTER devices are used for the anterograde placement of a voice prosthesis when replacing the voice prosthesis.

4. INDICATION

A voice prosthesis must not remain permanently in the shunt channel, but needs to be replaced at regular intervals, depending on its proper functioning.

The leakage of fluids with the transfer of saliva or food particles from the oesophagus to the trachea is a direct indication for an examination by the physician and, where necessary, replacement of the voice prosthesis used.

The **INSERTER** replacement is performed in patients who have been equipped with a voice valve following removal of the larynx (laryngectomy).

5. CONTRAINDICATION

The use of the **INSERTER** for voice prosthesis placement in patients with increased bleeding tendency, e.g. due to treatment with anticoagulants, must always be weighed carefully on account of the risk.

Special caution is required in patients following radiation therapy and/or chemotherapy, as there is an increased risk of puncture-induced complications due to possible tissue damage (e.g. scar tissue, fibroses and the like).

6. COMPLICATIONS AND ADVERSE REACTIONS

WARNING!

Should the voice prosthesis come loose from the shunt channel or during insertion of the placement instrument, there is a risk of aspiration.

In the event of aspiration, it must be attempted to expel the voice prosthesis through strong coughing.

If the voice prosthesis cannot be coughed out, it must immediately be removed surgically.

CAUTION!

When the voice prosthesis is inserted, this can lead to irritations, coughing or slight bleeding.

Further complications:

- General discomfort due to manipulation on the shunt channel
- Expansion/dilatation of the shunt channel
- Inflammatory reactions and formation of granulations in the region of the shunt channel

PLEASE NOTE!

Inappropriate handling of the **INSERTER** can cause material damage to the voice prosthesis.

7. PRODUCT DESCRIPTION

The **INSERTER** consists of a placement instrument (s. Fig. A) for the secure fixation of the voice prosthesis for insertion with the applicator (s. Fig. B) into the shunt channel.

There are two marking lines (s. Fig. E) on the applicator.

The applicator is mounted on an adapter (s. Fig. C).

NOTE:

In addition, also heed the instructions for use of the voice prosthesis used.

8. APPLICATION

1. Prior to removing the voice prosthesis to be replaced, check its seating in the shunt channel.
The correct length of the voice prosthesis is of essential importance for avoiding complications such as peri-prosthetic leakage, pressure necroses or the like.
2. Check the packaging for damages and check the expiry date.
Do not use the product after expiry of this date, or if the packaging and/or the product exhibit visible damages.
3. Hold the placement instrument such that the long slit is visible towards the distal end (s. Fig. 1).
4. Feed the safety thread of the Uniek voice prosthesis from the slotted side of the inserter halfway through the hole provided until the safety thread protrudes by 8-10 mm (s. Fig. 1).
5. Fixate the safety thread of the voice prosthesis through the slit by pulling slightly towards the distal end until it clicks noticeably while holding the loose end on the back (s. Fig. 2).
6. Check to ensure that the voice prosthesis is securely seated on the placement instrument.
The voice prosthesis is correctly secured when the distance between the top of the placement instrument and the tracheal flange is approx. 5 mm.
7. Pull on the voice prosthesis to check if it is securely in place.
8. For better lubrication of the applicator to advance the voice prosthesis, use a water-based lubricant by putting an approximately pea-sized amount into the funnel of the applicator.
9. Place the applicator in the adapter with the bevelled tip facing downwards (s. Fig. 3).
10. Push the oesophageal flange with the placement instrument forward into the applicator by rotating it slightly until the voice prosthesis touches the adapter (s. Fig. 4).
11. Then remove the loaded applicator together with the placement instrument from the adapter.
12. Make sure that the blue marking ring of the placement instrument (see Fig. D) is clearly visible between the two marking lines of the applicator (s. Fig. 5).
13. The voice prosthesis can now be introduced into the shunt channel using the applicator.
14. Gently insert the tip of the applicator into the shunt channel with a slight twisting motion (s. Fig. 6).
15. Slowly push the voice prosthesis beyond the 2nd marking line with the placement instrument (s. Fig. 7).
16. The oesophageal flange will open in the lumen of the oesophagus.
17. Now cautiously retract the applicator via the placement instrument (s. Fig. 8).
18. Pull gently on the placement instrument.
19. The tracheal flange opens in the lumen of the trachea.
20. The correct position of the voice prosthesis is verified by rotating the voice prosthesis with the placement instrument and by pulling slightly (s. Fig. 9).

21. When the tracheal flange lies in the fistula, it can be manipulated into the trachea using atraumatic forceps (s. Fig. 10).
22. Perform a phonation test and a swallowing test with the patient in order to verify the leakproofness.
23. If the voice prosthesis is in the correct position, cut the safety thread flush with the tracheal flange (s. Fig. 11).

9. HYGIENE INSTRUCTIONS

A **LARYVOX® INSERTER UNIEK PRO** must not be cleaned or disinfected. Sterilisation and reuse can impair the function and is therefore not permissible either.

10. STORAGE

Store in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

LARYVOX® INSERTER UNIEK PRO is a single-patient product and is intended for one single use only.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.


Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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LARYVOX® är ett i Tyskland och EU inregistrerat varumärke som ägs av Andreas Fahl Medizintechnik-Vertrieb GmbH, Köln.