
FAHL

# TRACHEOPORT® JUNIOR

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# CONTENT

<b>1.0</b>	<b>Introduction</b>	<b>26</b>	<b>7.0</b>	<b>Maintenance / reprocessing for reuse in another patient</b>	<b>37</b>
1.1	Notes on the instructions for use	26	7.1	Basic instructions	37
1.2	Function	27	7.2	Reprocessing for reuse	37
1.3	Purpose	27	7.2.1	Bacterial filter	38
1.4	Scope of delivery	28	7.2.2	Suction tube	38
1.5	Transport and storage	28	7.2.3	Secretion container	38
1.6	Legend of the pictograms	28	7.2.4	Lid of container	38
			7.2.5	Device surface	38
<b>2.0</b>	<b>Safety information</b>	<b>29</b>	<b>8.0</b>	<b>Troubleshooting</b>	<b>39</b>
<b>3.0</b>	<b>Set-up and start-up</b>	<b>31</b>	<b>9.0</b>	<b>Accessories, consumables</b>	<b>40</b>
3.1	Operating elements	31	9.1	Accessories	40
3.2	Connection to the power supply	31	9.2	Consumables	40
3.3	Start-up	31	<b>10.0</b>	<b>Technical data</b>	<b>41</b>
3.4	Charging the battery	31	<b>11.0</b>	<b>Disposal</b>	<b>42</b>
3.4.1	Charging with the power supply	31	<b>12.0</b>	<b>Instructions with regard to EMC</b>	<b>43</b>
3.4.2	Charging with the car connection cable	31	12.1	Guidelines and manufacturer's declaration - electromagnetic immunity	43
			12.2	Guidelines and manufacturer's declaration - electromagnetic immunity for TRACHEOPORT® JUNIOR	
<b>4.0</b>	<b>General operation</b>	<b>32</b>	12.3	Guidelines and manufacturer's declaration - electromagnetic emissions	45
4.1	Suction tube	32	12.4	Recommended safety distances between portable and mobile HF telecommunication devices and TRACHEOPORT® JUNIOR	45
4.2	Setting the vacuum	32	<b>13.0</b>	<b>Legal notices</b>	<b>46</b>
4.3	Suctioning	32			
4.4	Rinsing the tube	32			
<b>5.0</b>	<b>Operation</b>	<b>33</b>			
5.1	Secretion container and bacterial filter	33			
5.2	Inserting the bacterial filter and the secretion container	33			
5.3	Connecting the suction tube	33			
<b>6.0</b>	<b>Cleaning, disinfection, reprocessing</b>	<b>34</b>			
6.1	Basic instructions	34			
6.1.1	Bacterial filter	34			
6.1.2	Air filter	34			
6.1.3	Suction tube	34			
6.1.4	Secretion container	35			
6.1.5	Lid of container	35			
6.1.6	Device surface	35			
6.2	Oversuctioning	35			
6.3	Cleaning instructions	35			
6.4	Disinfection	36			
6.5	Recommended instrument disinfectants	36			
6.6	Recommended surface disinfectants	36			

To request further information, accessories, consumables or spare parts, contact:

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EN

## 1.1 NOTES ON THE INSTRUCTIONS FOR USE

These instructions for use contain important information as to how you can operate TRACHEOPORT® JUNIOR safely, properly and effectively. It helps avoid hazards and reduce repair costs and outage times, thus enhancing, amongst other features, the reliability and prolonging the service life of the device. They not only serve for briefing persons who operate the device, but are also intended to be used as reference.

**The instructions for use must always be kept available close to the device.**

Technically correct application ensures the operational safety and the functionality of TRACHEOPORT® JUNIOR and is therefore indispensable, besides cleaning performed on a regular basis.

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- The product TRACHEOPORT® JUNIOR bears the CE marking in accordance with EC Directive 93/42/EEC for medical devices and meets the basic requirements of Annex I of this directive.
- The quality management system applied at Andreas Fahl Medizintechnik-Vertrieb GmbH is certified according to the international standards DIN EN ISO 13485 .
- Read Section 2.0 „Safety information“ before taking the device into operation for the first time in order to avoid possible hazardous situations.
- These instructions for use correspond to the current state of the technical safety standards at the time of going to print.
- Reprinting – including in excerpts – is subject to the written consent of Andreas Fahl Medizintechnik-Vertrieb GmbH.
- Subject to alterations. Errors excepted.


**These instructions for use are valid for the following device:**

**TRACHEOPORT® JUNIOR (REF 63800)**

in connection with the accessory parts specified in Section 9.0, as well as for all devices of identical construction with special voltage.

**Please keep this document for later use!**

# 1.0 INTRODUCTION

## TRACHEOPORT® JUNIOR

- ❶ Vacuum gauge
- ❷ Vacuum control
- ❸ Battery capacity status indicator
- ❹ Operating status indicator
- ❺ Bacterial filter
- ❻ Secretion container
- ❼ Connecting tube
- ❽ Secretion outlet with stopper
- ❾ Connection for suction tube with stopper
- ❿ ON / OFF button
- ⓫ Connection for power supply unit

## 1.2 FUNCTION

- TRACHEOPORT® JUNIOR is a particularly handy small-size suctioning device.
- The device is powered by a maintenance-free electric swing-piston pump. When in operation, the pump generates a vacuum in the tubing system and in the secretion container, with the aid of which secretions or rinsing fluids are suctioned off (via a sterile suction catheter).
- The fluid is collected in the secretion container.
- The final vacuum, and thus also the suction power, can be adjusted with the aid of the fine control and the vacuum gauge.
- A bacterial filter prevents intrusion of bacteria into the inside of the unit.
- In addition, a mechanical overflow protector (float) integrated on the lid of the secretion container prevents accidental aspiration of secretions into the pump head.

## 1.3 PURPOSE

**Name:**  
TRACHEOPORT® JUNIOR

**Main function:**  
Temporary and spontaneous aspiration of secretions, blood and body fluids that typically arise during suctioning of the airways.

**Medical indication / application:**  
Suctioning of the upper airways

**Specification of the main functions:**  
Drainage and temporary collection of body secretions. A vacuum is generated by means of an electric suctioning pump. The integrated secretion container enables temporary collection of the body secretions that have been drained off.

**Application organ:**  
Upper airways (oral cavity, nasopharyngeal space and bronchial system)

**Duration of application:**  
Temporary application in patients undergoing suctioning of the airways (operating time of the battery up to 60 minutes)

**Application environment:**  
The application environment for airway suctioning is the clinical, medical practice, nursing and homecare field. Application must only be carried out by medically trained and briefed persons.

**Contraindication:**  
Not suitable for:

- continuous operation in connection with drainage in the low-vacuum range (e.g. thoracic drainage or wound drainage),
- long-lasting use in connection with endoscopy,
- use outside of the medical field,
- aspiration of flammable, caustic or explosive substances,
- suctioning in areas exposed to explosion hazards.

**TRACHEOPORT® JUNIOR is an active medical device.**

**Reprocessing:**  
The product is intended for repeated use. The accessories are partly reusable. For information on reprocessing, cleaning and disinfection, refer to instructions for use (see Section 6.0).

**TRACHEOPORT® JUNIOR must not be used:**

- for continuous operation in connection with drainage in the low-vacuum range (e.g. thoracic drainage) or long-lasting use in connection with endoscopy;
- outside of the medical field (for aspiration of flammable, caustic or explosive liquids/gases)







# 1.0 INTRODUCTION

EN

## 1.4 SCOPE OF DELIVERY

- This device was subjected to a thorough function check and was carefully packaged before being shipped. Please, nevertheless check the content of the packaging for completeness by comparing it with the scope of delivery immediately upon receipt. Besides the basic unit, the scope of delivery comprises the following:

### TRACHEOPORT® JUNIOR

	Secretion container with lid
	Bacterial filter
	Connecting tube for bacterial filter
	Air filters
	Suction tube, approx. 1.80 m long, incl. clip
	Fingertip
	Power supply / 2-pin 100-240 V~ 50-60 HZ
	Mains connection cable 2-pin, 230 V
	12 V car adapter
	Carrying bag
	Rinsing tank; 0.5 ml

## 1.5 TRANSPORT AND STORAGE

- The device must only be transported in the enclosed transport bag.
- Please promptly document and report transport damages.
- After being transported at temperatures below the freezing point, the device must be left to stand at room temperature

for up to six hours. If the device is not acclimatised, it must **not** be operated, as the membranes of the unit could be damaged.

### Ambient conditions:

Transport / storage: -20 to +60°C;  
10 to 95 % air humidity  
without condensation at an air  
pressure of 700 to 1060 hPa

Operation and charging  
of the battery:

+10 to +40°C;  
30 to 75 % air humidity  
without condensation at an air  
pressure of 700 to 1060 hPa

## 1.6 LEGEND OF THE PICTOGRAMS

Observe instructions for use!

Device of protection class II

Type BF applied part

Serial number

Catalogue number

Manufacturer

Please observe the laws and regulations  
which apply to the disposal of electrical  
and electronic devices in the respective  
countries.

Battery capacity status indicator

Operating status indicator

Not for reuse

Attention, important information

IPN<sub>1</sub>N<sub>2</sub>

IP classification

### GENERAL SAFETY INSTRUCTIONS

Please read these instructions for use through carefully before taking TRACHEOPORT® JUNIOR into operation.

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
- no original accessory parts and spare parts are used,
- the instructions regarding proper use in these instructions for use are disregarded,
- assembly, readjustments, modifications, extensions and repairs have not been carried out by persons and companies who have been authorised by Andreas Fahl Medizintechnik-Vertrieb GmbH.

The TRACHEOPORT® JUNIOR complies with the recognised rules of technology and the regulations of the EU Directive for Medical Devices 93/42 EEC.

No warranty claims are accepted in the event of damages or malfunctions that have arisen due to the use of accessories or consumables from external companies.

Only transparent tubes must be used.

The repeated use of components labelled with  is prohibited. If used repeatedly, these components use their function and there is a high risk of infection.

The device must not be operated in environments exposed to explosion hazard and enriched with oxygen.

### RISK OF INJURY!

TRACHEOPORT® JUNIOR is designed for aspirating body fluids, secretions and rinsing fluids that typically arise during suctioning of the airways. The device must not be used to aspirate any explosive, flammable or caustic gases or liquids.

**Do not use parts packaged in sterile packaging anymore if the packaging has been damaged during transport or storage → risk of infection for the patient!**

The device must not be operated in areas exposed to splash water or in zones exposed to explosion hazard.

If the device is to be fully ready for use at any time for emergencies, then it must be continuously charged by means of the battery charger while being stored.

Only use the power supply unit included in delivery.

This suctioning device must not be operated without disposable bacterial filter.

TRACHEOPORT® JUNIOR must only be used by persons who have been briefed in the medical application.

The suction tube to be used must never come directly into contact with the place where suction is to be applied. A sterile suction catheter or a medically approved suctioning set must always be used for suctioning.

Pull the plug of the power supply out of the wall socket to disconnect the device from the mains power supply in each of the following situations:

- before cleaning and performing maintenance on the device,
- before emptying the container,
- each time before leaving the room.

Do not pull the plug out of the wall socket by pulling on the cable! Never touch power supply/ mains connection cable with wet hands.

#### **Caution!**

There is a risk of strangulation due to the long suction tube, the mains plug of the power supply and the low voltage cable of the power supply unit.

## 2.0 SAFETY INFORMATION

EN

### EQUIPMENT DAMAGE!

The device must not be taken into operation:

- if the supply cable or the plugs are damaged,
- if the device has fallen down once,
- if the device shows obvious safety defects.

Clean the device and return it for repair via the dealer from whom you have purchased TRACHEOPORT® JUNIOR.

The device must be checked at regular intervals for proper function and technical safety defects, e.g. plug contacts, secretion container, housing, etc.

The ambient conditions specified in Section 1.5 'Transport and storage' must be observed.

Never immerse the device in water, even when the device is in non-activated state.

Before taking the device into operation, the device itself, the secretion container, accessories, connection cables and tubes must be checked for damages. Damaged cables and tubes must be replaced immediately. The device must be checked for proper function before use.

No liquid must be allowed to penetrate into the device. If liquid has penetrated into the device, then the device must only be taken into operation again after it has been checked by the customer service.

The device must be set up on a firm, level surface. If the device is set up on an uneven surface (e.g. mattress, cushion, upholstered chair), there is a risk that the device may tip.

Before the device is connected, it must be checked whether the mains voltage and frequency indicated on the device are the same as the voltage and frequency of the supply network.

Only use proper, regulation-compliant mains connections and extension cables. Avoid moisture on plug and switch unit.

### HANDLING OF RECHARGEABLE BATTERIES

The rechargeable battery must be fully charged before being taken into operation for the first time!

Deep discharge destroys the rechargeable batteries. **Therefore, charge the batteries every 3 months even when TRACHEOPORT® JUNIOR is not being used.**

Always store devices that are operated with rechargeable batteries in charged condition.

If the device has not been in operation for a longer period of time, the full capacity of the rechargeable batteries is only once more available after approx. 4 full charging and discharging cycles.

Have worn-out rechargeable batteries replaced immediately by Andreas Fahl Medizintechnik-Vertrieb GmbH or by authorised companies / persons. Mains operation of the device with worn-out rechargeable batteries can destroy the electronic charging system, or excessive power consumption can lead to spontaneous shut-offs of the device.

Heat destroys the rechargeable batteries. Therefore, avoid exposure to direct sunlight as well as the proximity of heaters. The German Federal Environment Agency (Umweltbundesamt) recommends an ideal ambient temperature between 10-25 °C. Ambient temperatures below -20 °C and above 50 °C decrease the lifespan of the battery.

If the available capacity (operating time) of the rechargeable battery has dropped to less than 80 % of the capacity of a new battery, the battery should be replaced by Andreas Fahl Medizintechnik-Vertrieb GmbH or by authorised companies / persons.

**Rechargeable batteries are wear parts. Due to the way they are manufactured, a maximum of 800-900 charging cycles can be achieved with the batteries, provided that they are stored, taken into operation and used properly, and in particular in accordance with the specifications contained in these instructions for use. If the batteries are not stored, handled and/or taken into operation properly, a service life considerably below 800-900 charging cycles can result.**

## 3.0 SET-UP AND START-UP

### 3.1 OPERATING ELEMENTS

- ❶ Vacuum gauge
- ❷ Vacuum control
- ❸ Battery capacity status indicator
- ❹ Operating status indicator
- ❺ Bacterial filter
- ❻ ON / OFF button
- ❼ Low voltage connection socket
- ❽ Low voltage cable of the power supply
- ❾ Mains plug of the power supply


Figure 2a


Figure 2b


Figure 3a


Figure 3b

Figure 4

### 3.2 CONNECTION TO THE POWER SUPPLY

Check whether the mains voltage and frequency indicated on the device are the same as the voltage and frequency of the supply network.

Check the connection cables for possible damage. If the cable shows damage, replace it immediately!

### 3.3 START-UP

- TRACHEOPORT® JUNIOR is delivered ready for operation.
- Lift the device out of the packaging. Check whether the voltage values indicated on the power rating plate agree with the voltage supplied.
- Always place the device on a level, secure surface.
- The rechargeable battery must be fully charged prior to first use. To this end, charge the device as described in Section 3.4!
- Switch the device on with the aid of the ON/OFF button (Figure 2b, ❻). If the operating status indicator (Figure 2a, ❹) lights up in green, the device is ready for operation. If the operating status indicator (Figure 2a, ❹) lights up in yellow, the rechargeable battery must be charged. (see Section 3.4).
- Be sure to heed the safety instructions in Section 2.0 before taking the device into operation for the first time.
- After being transported at temperatures below the freezing point, the device must be left to stand at room temperature for up to six hours. If the device is **not** acclimatised, it must **not** be operated, as the membranes of the unit could be damaged.
- Always have at least one bacterial filter ready, since it is not allowed to operate the device without this!

### 3.4 CHARGING THE BATTERY

The battery is charged via the 12 V low-voltage connection of the device. While charging is in progress, the battery capacity status indicator (Figure 2a, ❸) flashes green. As soon as the battery has been fully charged, the battery capacity status indicator (Figure 2a, ❸) lights in green. The device can be charged or operated with either the power supply (Figure 3a) or the car connection cable (Figure 3b).

#### 3.4.1 OPERATION / CHARGING WITH BATTERY CHARGER

Connect the low voltage cable of the power supply (Figure 3a, ❽) to the low voltage connection of the device (Figure 2b, ❷). Plug the mains plug of the power supply (Figure 3a, ❾) into the wall socket.

#### 3.4.2 OPERATION / CHARGING WITH CAR CONNECTION CABLE

Connect the low voltage connection of the device (Figure 2b, ❷) to the cigarette lighter socket of the car, using the special car cable (Figure 3b).

#### CAUTION!

When the power supply or the car connection cable is connected, the full suction power is available, even if the battery is empty. If the suctioning device is switched on via the ON/OFF button, the battery capacity status indicator goes out and the operating status indicator lights up in green.

## 4.0 GENERAL OPERATION

EN

Please read the preceding section before going on to read these instructions for use!

### IMPORTANT SAFETY INSTRUCTIONS!

**Caution: Do not perform suctioning in the airways until you have been briefed by specialist staff.**

Only use special sterile single-use suction catheters for suctioning. Always use a new suction catheter for each suctioning process.

Take care to ensure that the secretion container is emptied in good time. As soon as the secretion container is half full, the container must be emptied (this applies to all fields of application)!

If the secretion container is overfilled (which must be avoided), the overflow protector responds and the device does not suction anymore. Empty the container and replace the filter!

Check the vacuum gauge regularly!

If secretions have been aspirated into the pump due to incorrect use or manipulation, the device must be repaired by the manufacturer or by an authorised specialist company.

Use suitable suction catheters, suctioning extensions or suctioning sets for suctioning.

Keep an eye on the fluid level in the secretion container during the suctioning process.

**AUXILIARY AIR OPENING ④ OPEN =  
Suctioning process is interrupted  
(e.g. when inserting the suction catheter)**

**CLOSE AUXILIARY AIR OPENING ④  
WITH FINGER =  
Suction is applied**

### 4.3 SUCTIONING

- Now, insert the suction catheter as you have been shown by the specialist staff and start the suctioning process.
- Control the suctioning process by means of the auxiliary air opening (Figure 5, ④) on the Fingertip.
- The oversuctioning protector reliably prevents the penetration of fluid into the pump. You should nevertheless empty the container when it is half full though.

### 4.4 RINSING THE TUBE

- After each suctioning process, dispose of the suction catheter and briefly flush the suction tube with clean water or disinfectant (see Section 6). You will find it helpful to use a rinsing bottle in which clean water can be kept available.

### 4.1 SUCTION TUBE

- Connect the suction catheter (Figure 5, ③) to the suction tube (Figure 5, ①) with the aid of the Fingertip (Figure 5, ②).


Figure 5

- ① Suction tube
- ② Fingertip
- ③ Suction catheter
- ④ Auxiliary air opening

### 4.2 SETTING THE VACUUM

- Set your desired vacuum by keeping the auxiliary air opening of the Fingertip (Figure 5, ④) closed and allowing the vacuum to build up. Then turn up the vacuum control (Figure 2a, ②) until the vacuum gauge indicates the desired vacuum.
- Choose a suction catheter (Figure 5, ③) with the correct size or a suctioning set, which can be purchased via specialist dealers.
- Connect the suction catheter to the suction tube as described in Section 4.1.

## IMPORTANT SAFETY INSTRUCTIONS REGARDING THE CONTAINER SYSTEM

**The device must not be operated without the bacterial filter, oversuctioning protector and air filter! Therefore, have at least one spare bacterial filter and one spare air filter ready!**

Check each time before using the device whether the bacterial filter / oversuctioning protector / air filter and the container lid are clean and dry. Moist or contaminated filters must be replaced by new filters.

For your own safety, use disposable gloves when replacing the bacterial filter / air filter!

The bacterial filter / air filter is a single-patient product. The bacterial filter / air filter must therefore be replaced each time when the device is used for another patient. If the device is used exclusively for one patient, the filters must be replaced at least every 4 weeks (depending on frequency of use).

### 5.1 SECRETION CONTAINER AND LID OF CONTAINER

- Fit the container lid onto the secretion container, holding it horizontally while the container is standing on a firm surface.
- Press the container lid onto the container (Figure 6 / a) with both hands as far as it will go, exerting light pressure & close the lid with a clockwise rotation.
- Then insert the secretion container into the device (Figure 6 / b,c).
- Put the head on the secretion container (Figure 6 / d).

### 5.3 CONNECTING THE SUCTION TUBE

- Connect the suction tube to the connecting piece at the side of the secretion container lid (Figure 6 / f ③), turning it slightly.
- Likewise turn the tube slightly to remove it.
- Open the secretion container by turning it counterclockwise.

#### CAUTION!

**Failure to observe this sequence can lead to a reduction of suctioning power!**

### 5.2 CONNECTING THE BACTERIAL FILTER TO THE SECRETION CONTAINER

- The filter is connected between the connection on the suctioning device (Figure 6 / e ①) and the connection for the suctioning device in the centre of the secretion container lid (Figure 6 / e ②) with the aid of the two short connecting tubes.


Figure 6

## 6.0 CLEANING, DISINFECTION, REPROCESSING

Careful handling of the device has a decisive impact on the reliability and safety of the device. A daily visual inspection of the device for external damage and contamination is indispensable.

The cleaning measures described in the following must be carried out even when the device is not being used regularly. Regular replacement of the filters for TRACHEOPORT® JUNIOR at 4-week intervals and function checks at intervals of not more than 3 months are mandatory in order to ensure the functionality of the device.

These hygiene measures are necessary measures to protect patients and users and to maintain the functional reliability of the suctioning device.

If the device is to be used for another patient or if oversuctioning of the device has occurred, professional reprocessing at the manufacturer's establishment or at an authorised specialist company in accordance with the Medical Devices Operator Ordinance (MPBetreibV), the Medical Devices Act (MPG) and the Federal Association for Medical Technology (BV-Med) guideline is required.

### CAUTION!

**The suctioning device is a single patient product and must only be used for one patient, since otherwise there is a risk of infection. This suctioning device can be reprocessed for reuse in another patient. When doing this, observe the instructions in this regard provided in the instructions for use, particularly in Section 7. (Maintenance / reprocessing for reuse in other patient). The manufacturer accepts no liability for unauthorised and incorrect reprocessing of this device by external companies!**

### 6.1 BASIC INSTRUCTIONS

**The suctioning device must be cleaned in accordance with the following cleaning and maintenance instructions.**

The following instructions are the result of long years of experience. Shorter cycles may be required, depending on the application and the users' experience.

- A daily visual inspection of the device for external damage and contamination is indispensable.
- We recommend that you document all maintenance and replacement processes in writing as a matter of principle.
- Depending on the field of application of the device, e.g. in the nursing field, additional precautions must be taken before starting with suctioning, in order to ensure that the safety of the patient, the user and third parties is not jeopardised. These include for instance thorough disinfection of hands and/or wearing of disposable gloves.
- Disconnect the mains connection cable from the device before cleaning the device.
- The cleaning measures described do not replace the respective regulations valid for the establishment!
- Observe the instructions for use of the manufacturers of the disinfectants, particularly the information regarding concentration and material compatibility.
- Disinfection of the device should only be performed with the agents specified by the manufacturer, in order to avoid damage to the device.
- **CAUTION!**  
**Some disinfectant solutions stain the lid components and connecting tubes somewhat. This has no impact on the properties of the materials though.**
- Basically, the reprocessing instructions given in Section 7 must be observed when the device is used for another patient.
- Partial cleaning in an automatic washer-disinfector (secretion container and container lid) is in principle also possible. Thermal disinfection is performed at a temperature of 93°C.

- All product parts that come into contact with bodily fluids (secretion container and cap, suction tube – with or without Fingertip, rinsing tank) can also be thermally disinfected in the domestic environment by means of boiling.

#### **Instructions for thermal disinfection in domestic environments:**

1. Place a sufficiently large saucepan on the hob and fill with water.
2. Bring the water to the boil.
3. Then add the inhalation mask to the boiling water.  
Make sure that the objects to be disinfected are completely covered with water and that all surfaces are rinsed with water.  
Stir the water occasionally with a wooden spoon or similar.
4. Boil the items for 10 minutes in continuous boiling water.
5. Then remove the objects from the water using a wooden spoon, grill tongs or similar  
**Caution! Danger of scalding!**
6. Let the objects cool down and then dry them with a clean, lint-free cloth.

### CAUTION!

**Boiling may cause discolouration of the material. However, this has no effect on the functionality of the inhalation masks.**

- Care must be taken to ensure that the suctioning device and the suction tubes are stored under hygienic conditions!

### 6.1.1 BACTERIAL FILTER

The bacterial filter prevents microorganisms and secretions from penetrating into the device and from being blown out of the device, and thus serves to protect the users and the device.

- The filter cannot be cleaned or disinfected
- For reasons of hygiene, we recommend replacement at least every 4 weeks. If the filter is strongly contaminated, it must be replaced immediately.
- In order to prolong the service life of the bacterial filter, it is advisable to fill the secretion container only to approximately half its capacity. Always use an original bacterial filter.
- **CAUTION!**  
**The suctioning device must not be operated without bacterial filter.**
- Make sure you always have a sufficient number of spare filters available.

### 6.1.2 AIR FILTER

The air filter (see page 15 figure 9) prevents microorganisms and secretions from being blown out of the device, and thus serves to protect the users and the device.

- The filter cannot be cleaned or disinfected
- For reasons of hygiene, we recommend replacement at least every 4 weeks. If the filter is strongly contaminated, it must be replaced immediately. If the air filter is moist, or if you detect deposits, please observe Section 6.2 of these instructions.
- In order to prolong the service life of the air filter, it is advisable to fill the secretion container only to approximately half its capacity. Always use an original air filter.
- **CAUTION!**  
**The suctioning device must not be operated without air filter.**
- Make sure you always have a sufficient number of spare filters available.

# 6.0 CLEANING, DISINFECTION, REPROCESSING

## 6.1.3 SUCTION TUBE

- In order to prevent encrustation of the suction tube with dried secretions, the tube must be flushed with clear water after each suctioning process. The water can be sucked through the fingertip into the enclosed rinsing tank with the help of the suction device. Please only fill the rinsing tank to about half its capacity.
- We additionally recommend disinfection with one of the instrument disinfectants recommended on page 36 at least once every day. Observe the respective instructions for use of the disinfectants!
- To dry the suction tube, it is advisable to hang it up. The clip included in the first delivery can be used for this.
- Frequent cleaning / disinfection / sterilisation can cause the suction tube to become discoloured and the material to become brittle. We therefore recommend replacement at least every 4 weeks.

## 6.1.4 SECRETION CONTAINER

- For reasons of hygiene, we recommend to empty the secretion container after each suctioning process and then to clean it vigorously with warm water, preferably with commercially available dishwashing detergent added. Remove tenacious contaminations with a commercially available bottle brush.
- **For instructions on how to remove the container, refer to Section 5.0 ,Operation‘ on page 33.**
- We recommend to disinfect the secretion container once every day with one of the instrument disinfectants recommended on page 36. Please note that the secretion container should first be emptied and rinsed with clear water in order to ensure a better cleaning result.

**The respective instructions for use of the disinfectants must be observed!**

**Note: Never hold the container by the lid to carry it. The filled container can come off the lid and fall down to the floor (Figure 7). Use the secretion outlet to empty the secretion container.**


Figure 7

## 6.1.5 LID OF CONTAINER

**The bacterial filter must ALWAYS be removed before cleaning.**

For instructions on how to remove the secretion container lid, refer to Section 5.0 on page 33.

- We recommend to disassemble the container lid after each suctioning process and to rinse the individual components (lid, float, tube connector) thoroughly under clear water. The lid must be absolutely dry before it is used again. When assembling the lid, care must be taken to ensure that the overflow protector is not damaged
- In order to achieve a better cleaning effect, first rinse the individual components for at least 10 seconds under clear water. Then disinfect the individual components with one of the instrument disinfectants recommended on page 36. **The respective instructions for use of the disinfectants must be observed!**

## 6.1.6 DEVICE SURFACE

- In the event of contamination, but at least once a week, the entire surface of the device should be cleaned with a moist cloth (never use a wet cloth) and should then be disinfected with one of the surface disinfectants listed on page 36.
- **Some disinfectants can cause discolouration of the material, which does not affect the functionality of the device though.**
- **Make sure that the device has first been disconnected from the mains power supply!**
- **The device must never be autoclaved, rinsed under running water or immersed in liquids!**

**If liquid has penetrated into the device, then the device must only be taken into operation again after it has been checked by the authorised customer service.**

## 6.2 OVERSUCTIONING

**When has a suctioning device been exposed to oversuctioning?**

A suctioning device is considered to have been exposed to oversuctioning if suctioning fluid has been suctioned into the inside of the device.

**How can I discern that oversuctioning has occurred?**

TRACHEOPORT® JUNIOR has an air filter on its underside (page 37 / Figure 9). To check, remove the cover cap. If moisture, gross contamination or deposits are visible in the air filter, this shows that the suctioning device has been exposed to oversuctioning. As a general rule, reduced suction power is a sign that the device has possibly been exposed to oversuctioning. If oversuctioning has occurred, the suctioning device must be reprocessed by the manufacturer or a certified partner of Andreas Fahl Medizintechnik-Vertrieb GmbH.

**A device that has been exposed to oversuctioning constitutes a hazard both for the patient and for the carer. For this reason, we recommend that the air filter should be checked on a regular basis.**

## 6.3 CLEANING INSTRUCTIONS

To improve the cleaning effect, commercially available dishwashing detergent can be added to the warm rinsing water. If secretion container and container lid are strongly contaminated, they should be left to soak in the rinsing water for a longer period of time. Tenacious contaminations must be removed with a soft brush or a soft cloth.

- **Partial cleaning in an automatic washer-disinfector (secretion container and container lid) is in principle also possible.**
- **Thermal disinfection is performed at a temperature of 93°C.**

The cleaning and maintenance instructions are the result of long years of experience. Shorter cycles may be required, depending on the application and the users' experiences.

**The following consumables must be replaced in accordance with the instructions given in Section 6.1 (Figure 8):**

- **The device must be disconnected from the mains for cleaning!**
- **The instructions specified in the instructions for use must be observed, particularly as regards the recommended agents.**

## 6.0 CLEANING, DISINFECTION, REPROCESSING

EN


Figure 8

### 6.4 DISINFECTION

The different device components can be disinfected with specific disinfectants:

- Suction tube (instrument disinfectant)
- Secretion container (instrument disinfectant)
- Container lid (instrument disinfectant)
- Device surface (surface disinfectant)

**In this regard, observe the instructions already mentioned in this section as well as the disinfectant manufacturers' instructions for use.**

Thorough cleaning should always be performed prior to disinfection, if such is required. As a general rule, disinfection is advisable to prevent cross-infection and is also advisable if the device is used in inpatient situations (e.g. hospitals, nursing homes and/or other health care facilities), in order to limit the risk of infections.

### 6.5 RECOMMENDED INSTRUMENT DISINFECTANTS

Disinfectant	Ingredients	(in 100 g)	Manufacturer
GIGASEPT FF neu (application concentrate)	Succinic acid dialdehyde Dimethoxy tetrahydrofurane Phosphonates, anionic tensides, non-ionic tensides, scents, methyl isothiazolinone	11.9 g 3.2 g	Schülke & Mayr GmbH, Norderstedt
Sekusept Plus (application concentrate)	Glucoprotamin	25.0 g	Ecolab Deutschland GmbH, Düsseldorf not for rinse container

### 6.6 RECOMMENDED SURFACE DISINFECTANTS

Disinfectant	Ingredients	(in 100 g)	Manufacturer
Dismozon pur steril (application solution)	Magnesium monoperoxyphthalate hexahydrate	80 g	BODE Chemie GmbH, Hamburg
Kohrsolin FF (application solution)	Glutaral Benzyl-C12-18-alkyldimethyl ammonium chlorides Didecyl dimethyl ammonium chloride	5 g 3 g 3 g	BODE Chemie GmbH, Hamburg
Mikrozid sensitive wipe	Quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl, chlorides Didecyl dimethyl ammonium chloride Quaternary ammonium compounds, benzyl-C12-14-alkyl [(ethyl-phenyl)methyl]dimethyl, chlorides	0.26 g 0.26 g 0.26 g	Schülke & Mayr GmbH, Norderstedt
Perform (application solution)	Pentapotassium bis(peroxymonosul- phate) bis(sulphate) Anionic tensides, non-ionic tensides, phosphonates, soap, scents	45 g	Schülke & Mayr GmbH, Norderstedt
Bacillol 30 Foam	Propane-1-ol Propane-2-ol Ethanol	45 g 25 g 4.7 g	BODE Chemie GmbH, Hamburg
Mikrobac forte	Benzyl-C12-18-alkyldimethyl ammonium chlorides N-(3-aminopropyl)-N-dodecylpropane- 1,3-diamine	19.9 g 5 g	BODE Chemie GmbH, Hamburg

If disinfectants containing aldehyde and disinfectants containing amines are used on the same object, this can cause discolourations.

# 7.0 MAINTENANCE / REPROCESSING FOR REUSE IN ANOTHER PATIENT

## 7.1 BASIC INSTRUCTIONS

For TRACHEOPORT® JUNIOR, there are no mandatory prescribed specifications in terms of times or time limits for maintenance as per Medical Devices Operator Ordinance beyond the requirements for cleaning/disinfection specified in the instructions for use. If reprocessing for reuse is performed, technical function checks as specified in the Medical Devices Operator Ordinance are required. In the commercial field, checks performed on a regular basis may be required in the scope of accident prevention regulations that apply at the workplace. For non-stationary electrical equipment, these checks are due at least every 24 months. Information in this regard can be obtained from the competent occupational insurance associations.

Regular thorough cleaning and disinfection of the tubes and the application parts, as well as operation of the device in accordance with the instructions for use, is assumed.

Check the air filter on the underside of the device on a regular basis. To do this, withdraw the cover and examine the colour of the air filter. If you discover discolourations/deposits/moisture: send the device to the manufacturer or to an authorised specialist company.

If the device is to be used for another patient, or if oversuctioning of the device has occurred, professional reprocessing at the manufacturer's establishment or at an authorised specialist company in accordance with the Medical Devices Operator Ordinance (MPBetreibV), the Medical Devices Act (MPG) and the Federal Association for Medical Technology (BV-Med) guideline is required.

For your partners and customers, Andreas Fahl Medizintechnik-Vertrieb GmbH offers unproblematic, fast reprocessing and checking of Fahl suctioning devices.

If the device is to be used for another patient or if it changes owner, the device must be reprocessed in accordance with the Medical Devices Operator Ordinance, the Medical Devices Act and the Federal Association for Medical Technology guideline. For a contaminated device that has been exposed to oversuctioning to be reused, it must first be repaired by the manufacturer, by a certified partner of Andreas Fahl Medizintechnik-Vertrieb GmbH or by a specialist dealer who has been authorised by Andreas Fahl Medizintechnik-Vertrieb GmbH.

For partners and customers, Andreas Fahl Medizintechnik-Vertrieb GmbH offers unproblematic, fast reprocessing and examination / technical safety checks for Fahl suctioning devices.

If the device is to be used for another patient, all parts that come into contact with suctioning fluid (container, lid, overflow protector) must, as a matter of principle, be cleaned and then disinfected. See in this regard the individual cleaning instructions (see Section 6.3). In addition, filters, suction tubes, Fingertip etc. must be replaced.

- Perform a visual inspection of the device, including tubes, secretion container and device connection cable, each time before you use the device.

**Replace damaged cables and tubes immediately!**

- Maintenance, or opening and repair of the device, must only be performed by the manufacturer or by an authorised specialist company. When doing this, technical and hygiene-related precautions, the safety instructions as well as the servicing instructions for TRACHEOPORT® JUNIOR must be observed.

- This device can be returned for repair to Andreas Fahl Medizintechnik-Vertrieb GmbH either directly or via the dealer from whom you have purchased the device.
- Before the device is sent in, all secretion container and tube components must be cleaned and then disinfected. The device itself must be subjected to a surface disinfection
- Andreas Fahl Medizintechnik-Vertrieb GmbH neither warrants faultless function, nor does Andreas Fahl Medizintechnik-Vertrieb GmbH accept liability for damage to persons or property, if
  - no original parts are used,
  - the instructions regarding proper use in these instructions for use are disregarded,
  - assembly, readjustments, modifications, extensions and repairs are carried out by specialists who have not been authorised by Andreas Fahl Medizintechnik-Vertrieb GmbH.
- No warranty claims are accepted in the event of damages or malfunctions that have arisen due to the use of accessories or consumables from external companies.
- The regulations and instructions that apply to the respective field of application must be observed.

## 7.2 REPROCESSING FOR REUSE

The handling of the device has a decisive impact on the reliability and safety of the device. The hygiene measures described in the preceding sections are necessary to protect patients and users and to maintain the functional reliability of the suctioning device. They do not replace reprocessing by the manufacturer or an authorised specialist.

**Before being passed on to a new patient/user, the suctioning device must be reprocessed in accordance with the manufacturer's specifications.**

### How can you identify a contaminated suctioning device?

The condition of the air filter is one (but not the only) possible sign pointing to contamination. Perform a visual inspection of the air filter (on the device underside) (Figure 9). To do this, you must remove the cover. If the air filter is moist or contaminated, or if it shows deposits, the device has been exposed to oversuctioning and is contaminated.

**In this case, the suctioning device must be repaired by the manufacturer or by an authorised specialist company. Special inspection and reprocessing measures are required for this.**


Figure 9

**Air filter for quick visual inspection for possible contamination**

### CAUTION!

**Even if you cannot ascertain that the device has been exposed to oversuctioning, the device must nevertheless be repaired by the manufacturer or by an authorised specialist company. The device can thereafter be used again.**

## 7.0 MAINTENANCE / REPROCESSING FOR REUSE IN ANOTHER PATIENT

### 7.2.1 BACTERIAL FILTER

- The bacterial filter must be replaced each time when the device is used for another patient.

### 7.2.2 AIR FILTER

- The air filter must be replaced each time when the device is used for another patient.

### 7.2.3 SUCTION TUBE AND FINGERTIP

- The suction tube and the tube connector (Fingertip) must be replaced each time when the device is used for another patient.

### 7.2.4 SECRETION CONTAINER

- Each time when the device is to be used for another patient, the secretion container must be disinfected with one of the instrument disinfectants recommended on page 36. Please note that the secretion container should first be emptied and rinsed with clear water in order to ensure a better cleaning result.

**The respective instructions for use of the disinfectants must be observed!**

### 7.2.5 LID OF CONTAINER

- Each time when the device is to be used for another patient, the lid of the container must be disinfected with one of the instrument disinfectants recommended on page 36. Please make sure that the container lid has first been disassembled into its individual components (lid, float and tube connector).

In order to achieve a better cleaning effect, first rinse the individual components (lid and float) for at least 10 seconds under clear water. Then disinfect the individual components.

**The respective instructions for use of the disinfectants must be observed!**

In order to prevent infections, always use disposable gloves during the entire procedure.

For instructions on how to remove the secretion container lid, refer to Section 5.0 on page 33.

### 7.2.6 DEVICE SURFACE

- Each time when the device is to be used for another patient, the entire surface of the device must be cleaned with a moist cloth (never use a wet cloth) and must then be disinfected with one of the surface disinfectants recommended on page 36.

Some disinfectants can cause discolouration of the material, which does not affect the functionality of the device though.

Make sure that the device has first been disconnected from the mains power supply!

The device must never be autoclaved, rinsed under running water or immersed in liquids!







## 8.0 TROUBLESHOOTING

This TRACHEOPORT® JUNIOR device has been subjected to a thorough inspection at the factory. Should a malfunction nevertheless arise, you may possibly be able to rectify it yourself if you observe the following instructions.

Symptom of malfunction	Possible cause	Remedy
<b>Device fails to start operating</b>	Rechargeable battery discharged	Connect the power supply. To enable the device to be operated without being connected to the mains power supply, the battery should be charged for 1-2 hours.
	Mains plug of the power supply is not well seated	Check all plug connections. Observe the indicator light; this must be lit if connection has been correctly established.
	Vacuum in the container	De-aerate container
<b>Insufficient power</b>	Rechargeable battery discharged	Charge battery
	Filter is obstructed	Replace filter
	Leak in the tubing or the secretion container lid	Check secretion container lid and tubing for firm seating
<b>1. Weak vacuum or no vacuum is indicated</b>  <b>2. Strong vacuum is indicated</b>	1.1 Bacterial filter is missing	Insert bacterial filter
	1.2 Leak in the tubing or on the secretion container lid	Check connections of connecting tube/ bacterial filter Check assembly position of the secretion container lid
	1.3 Secretions or blood have been aspirated into the device and the valve platelets of the unit are sticking	Check secretion container lid and tubing for firm seating If this is the case, the device must be sent in for repair
	2.1 Bacterial filter is obstructed	Replace bacterial filter
	2.2 The float of the overflow protector is blocking the secretion container lid inlet	Check inlet of secretion container lid. If necessary, empty container, clean overflow protector and check float for free movement

## 9.0 ACCESSORIES, CONSUMABLES

EN

	<b>Secretion container TRACHEOPORT® JUNIOR</b> 600 ml (1ml = 1cc)
	<b>Mains connection cable, 2-pin, 230 V for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR</b>
	<b>VE Electronic power sup- ply for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR</b>
	<b>Car connection cable for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR</b>
	<b>Air filters</b>
	<b>Bacterial filter for suctioning devices</b>
	<b>TRACHFLOW® LINE PRO</b> approx. 1.30 m long
	<b>Rinsing tank; 0.5 ml</b>

9.1 ACCESSORIES	REF
TRACHEOPORT® JUNIOR / secretion container, 600 ml	63800-1
Mains connection cable, 2-pin, 230 V for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR	60602
VE Electronic power supply for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR	60603
Car connection cable for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR	60604

9.2 CONSUMABLES	REF
Bacterial filter for TRACHEOPORT® COMPACT/ TRACHEOFIRST® COMPACT/ TRACHEOPORT® JUNIOR	60800-01
Air filter, pack of 10	50000-07
TRACHFLOW® LINE PRO suction tube with integrated Fingertip approx. 1.30 m long	60506
<b>Stoma catheters, pack of 30</b> sterile, 50 cm long, 06 Ch – 18 Ch, straight, central aperture, 2 small lateral aper- tures, packed individually in sterile packaging	<b>68000</b>
68000-06	green 06 Ch
68000-08	blue 08 Ch
68000-10	black 10 Ch
68000-12	white 12 Ch
68000-14	green 14 Ch
68000-16	orange 16 Ch
68000-18	red 18 Ch

# 10.0 TECHNICAL DATA

Suction capacity of unit	24 ± 2 l/min
max. vacuum at sea level	-76 kPa* (-760 mbar; -570 mmHg) ± 4 kPa
Vacuum indicator	-1...0 bar (± 25 mbar) (mm Hg; kPa)*
Vacuum control	continuously adjustable -20 kPa* (min.) to -75 kPa* (max.)
Auxiliary air control	manually via Fingertip
Secretion container	600 ml secretion container
Suction tube	ø 6 mm, approx. 1.80 m long
Power supply (power supply)	100-240 Vac (+/-10%) 50/60 Hz
Low voltage power supply	12 V DC ± 10 %
Operating time with battery (without mains power supply)	Operating time of battery approx. 60 min Permanent operation
Operating time with 12V DC mains power supply (via power supply or 12V car supply)	Permanent operation
Emergency operation	If the battery has been completely discharged, the device can operate with mains power supply
Current consumption (power supply)	up to 4 A
Power consumption	48 VA
Rechargeable battery	11.1 V; 4.4 Ah; Li ions; at least 500 charging cycles, with partial discharge correspondingly more; up to 800 - 900 charging cycles in total
Charging time	approx. 4 hrs.
Noise level	55.0 dB (A) @ 1 m
Ambient conditions Transport/storage	-20...+60°C 10...95 % air humidity without condensation at an air pressure of 700...1060 hPa
Operation	+10...+40°C 30...75 % air humidity without condensation at an air pressure of 700...1060 hPa
Dimensions, HxWxD	370 x 290 x 150 mm
Weight	2.4 kg incl. battery
Protection class (EN 60601-1)	II
Applied Part	Type BF
Protection type	IP 21

\* 1 bar = 750.06 mm Hg = 1000 hPa / depending on prevailing air pressure and ambient conditions.

Technical data last updated on: 2015-07-01

# 11.0 DISPOSAL

---

EN

- TRACHEOPORT® JUNIOR does not comprise any dangerous goods.
- The housing material is fully recyclable.
- Device and accessories must be decontaminated prior to disposal, since pathogen-containing secretion residues can lead to a hazard.
- Observe careful separation of materials.
- Observe country-specific disposal regulations (e.g. waste incineration).

## **Disposal in the EU**

The suctioning device described above is a high-quality medical device with a long service life. After the end of a life cycle, the device must be transferred for technically correct disposal. According to the EU directives (WEEE and RoHS), the device must not be disposed of via the general household waste. Please observe the existing laws and regulations for disposal of old equipment which apply in the respective country.

## **Disposal in Germany**

In the Federal Republic of Germany, the procedure for disposal of old electric equipment is specified in the Electrical and Electronic Equipment Act (ElektroG). The TRACHEOPORT® Junior is a type of product that is exempt from the German Electrical and Electronic Equipment Act (ElektroG) according to regulations of the German foundation for registering used appliances (Stiftung Elektro-Altgeräte Register). In order to ensure that the regulation-compliant disposal procedure is observed, either return the old device to your responsible specialist dealer, or send it directly for technically correct disposal to Andreas Fahl Medizintechnik-Vertrieb GmbH.

**Prior to disposal or transport, all accessory and tube components must be thoroughly cleaned and disinfected. The device itself must be subjected to a surface disinfection.**

# 12.0 INSTRUCTIONS WITH REGARD TO EMC

## 12.1 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

TRACHEOPORT® JUNIOR is intended for operation in the electromagnetic environment specified below. The customer or user of the TRACHEOPORT® JUNIOR device should ensure that the device is used in such an environment.

Electromagnetic immunity tests	IEC 60601- test level	Accordance level	Electromagnetic environment - guidelines
Conducted disturbances as per IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz	10 V	Portable and mobile radio devices should not be used at a distance closer to the TRACHEOPORT® JUNIOR device, including the cables, than the recommended safety distance as calculated with the equation that applies to the transmission frequency.  <b>Recommended safety distance:</b>  $d = 0.35 \sqrt{P}$  $d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz  $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz  where P is the nominal power of the transmitter in watt (W) according to the specifications of the manufacturer of the transmitter and d is the recommended safety distance in metres (m).  The field strength of stationary radio transmitters should be lower than the accordance level (b) for all frequencies as determined by an on-site investigation (a).  Interference is possible in the vicinity of devices marked with the following pictogram.  
Radiated HF disturbances as per IEC 61000-4-3	3 V / m 80 MHz to 2.5 GHz	10 V/m	

**NOTE 1**  
At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2**  
These guidelines may not be applicable in all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.

a The field strength of stationary transmitters such as base stations of wireless phones and land mobile radio units, amateur radio stations, AM and FM radio and television transmitters, cannot be predetermined exactly in theory. To establish the electromagnetic environment in regard to stationary transmitters, a study of the location should be considered. If the field strength measured at the location where the TRACHEOPORT® JUNIOR is used exceeds the above accordance levels, the TRACHEOPORT® JUNIOR should be observed in order to verify that the device functions as intended. If unusual performance characteristics are observed, additional measures may be required, such as for instance a change of orientation or a relocation of the TRACHEOPORT® JUNIOR device.

b Above the frequency range of 150 kHz to 80 MHz, the field strength should be below 3 V/m.

# 12.0 INSTRUCTIONS WITH REGARD TO EMC

EN

The TRACHEOPORT® JUNIOR device must not be used in an arrangement where it is placed directly next to other devices or stacked with other devices. If the device has to be operated close to or stacked with other devices, the device should be observed in order to verify its proper use when used in this arrangement

## 12.2 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY FOR TRACHEOPORT® JUNIOR

TRACHEOPORT® JUNIOR is intended for operation in the electromagnetic environment specified below. The customer or user of the TRACHEOPORT® JUNIOR device should ensure that the device is used in such an environment.

Electromagnetic immunity tests	IEC 60601- test level	Accordance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood or concrete or tiled with ceramic tiles. If the floor is covered with synthetic material, the relative air humidity must be at least 30 % betragen.
Fast transient electric disturbances/bursts as per IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input and output lines	± 2 kV for power supply lines	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment.
Surge voltages (surges) as per IEC 61000-4-5	± 1 kV common mode voltage ± 2 kV series mode voltage	± 1 kV common mode voltage	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment.
Voltage dips, short interruptions and variations of the power supply voltage as per IEC 61000-4-11	< 5 % $U_T$ (> 95 % dip of $U_T$ for 0.5 cycle)  40 % $U_T$ (60 % dip of $U_T$ for 5 cycles)  70 % $U_T$ (30 % dip of $U_T$ for 25 cycles)  < 5 % $U_T$ (> 95 % dip of $U_T$ for 5 s)	< 5 % $U_T$ (> 95 % dip of $U_T$ for 0.5 cycle)  40 % $U_T$ (60 % dip of $U_T$ for 5 cycles)  70 % $U_T$ (30 % dip of $U_T$ for 25 cycles)  < 5 % $U_T$ (> 95 % dip of $U_T$ for 5 s)	The quality of the supply voltage should be equivalent to that of a typical home or hospital environment. If the TRACHEOPORT® JUNIOR user requires continued functions even when the power supply is interrupted, it is recommended to power the TRACHEOPORT® JUNIOR from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	30 A/m	Magnetic fields at the mains power supply frequency should be equivalent to the typical values found in the business and hospital environment.

**NOTE**  $U_T$  is the alternating mains power supply voltage before application of the test levels.

## 12.0 INSTRUCTIONS WITH REGARD TO EMC

- Medical electrical devices are subject to special precautions with regard to EMC and must be installed in accordance with the EMC instructions specified below.
- Portable and mobile HF communication equipment can influence medical electrical devices.
- The use of accessories, transducers or cables other than those specified can lead to increased emission or a reduced electromagnetic immunity of the device or system.

### 12.3 GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

TRACHEOPORT® JUNIOR is intended for operation in an environment as that specified below. The customer or user of the TRACHEOPORT® JUNIOR device should ensure that the device is operated in such an environment.

Interference emission measurements	Accordance level	Electromagnetic environment - guideline
HF emissions as per CISPR 11	Group 1	TRACHEOPORT® JUNIOR uses HF energy exclusively for its internal function. Therefore, its HF emission is very low and it is unlikely that neighbouring electronic equipment will be disrupted.
HF emissions as per CISPR 11	Class B	TRACHEOPORT® JUNIOR is suitable for use in all facilities, including facilities in living quarters and facilities which are directly connected to a public power supply that also supplies buildings used for residential purposes
Emission of harmonics as per IEC 61000-3-2	Not applicable	
Emission of voltage fluctuations/ flicker as per IEC 61000-3-3	Not applicable	

### 12.4 RECOMMENDED SAFETY DISTANCES BETWEEN PORTABLE AND MOBILE HF TELECOMMUNICATION DEVICES AND TRACHEOPORT® JUNIOR

Recommended safety distances between portable and mobile HF telecommunication devices and TRACHEOPORT® JUNIOR			
TRACHEOPORT® JUNIOR is intended for operation in an electromagnetic environment with controlled HF disturbances. The customer or the user of the TRACHEOPORT® JUNIOR can help avoid electromagnetic interference by keeping the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the TRACHEOPORT® JUNIOR device – depending on the output power of the communication device as specified below.			
	Safety distance in metres depending on the transmission frequency		
Nominal power of the transmitter W	150 kHz to 80 MHz $d = [0.35] \sqrt{P}$	80 MHz to 800 MHz $d = [0.35] \sqrt{P}$	800 MHz to 2.5 GHz $d = [0.70] \sqrt{P}$
0.01	0.035	0.035	0.07
0.1	0.11	0.4	0.22
1	0.35	0.35	0.7
10	1.1	1.1	2.2
100	3.5	3.5	7.0
For transmitters of which the maximum nominal power is not indicated in the table above, the recommended safety distance d in metres (m) can be determined using the equation that belongs to the respective column, P being the maximum nominal power of the transmitter in watt (W) as specified by the manufacturer.			
<b>NOTE 1</b> At 80 MHz and 800 MHz, the higher frequency range applies.			
<b>NOTE 2</b> These guidelines may not be applicable in all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.			

## 13.0 LEGAL NOTICES

---

EN

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to TRACHEOPORT® JUNIOR, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to TRACHEOPORT® JUNIOR themselves caused thereby and to any consequential damages caused thereby.

If TRACHEOPORT® JUNIOR is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

TRACHEOPORT® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany.

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