





LARYVOX® PAD CONNECT

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I. FOREWORD

These instructions for use are valid for LARYVOX® PAD CONNECT.

The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE

LARYVOX® PAD CONNECT is used to align the skin level of an irregular tracheostoma that is situated below the skin level.

The speech option is improved by reducing undesired air loss.

LARYVOX® PAD CONNECT enables the connection of filter and valve systems via an additional 22 mm connector.

The product is intended for domestic use only.

III. WARNINGS

LARYVOX® PAD CONNECT is a single-patient product and is thus intended exclusively for one-time or repeated use in one individual patient.

When LARYVOX® PAD CONNECT is used for the first time, its suitability for sealing the tracheostoma should be checked.

In problematic individual cases, sealing of the tracheostoma may not be sufficient.

Should you experience skin irritation (e.g. redness, itchiness), discontinue the product immediately and consult a physician. This may be due to special skin sensitivity or allergy-related incompatibility.

IV. CONTRAINDICATIONS

Do not use LARYVOX® PAD CONNECT in case of allergy or intolerance to the material.

Do not wear LARYVOX® PAD CONNECT during laser treatment or electrocautery.

There is a considerable risk of injury from exposure to heat and, if applicable, from the release of toxic combustion products.

Use of the LARYVOX® PAD CONNECT should be avoided while undergoing radiotherapy due to the risk of skin irritations.

After undergoing radiotherapy, only use LARYVOX® PAD CONNECT after consulting your treating physician. Start using the tape at the earliest 2 weeks after completing radiotherapy.

V. PRODUCT DESCRIPTION

LARYVOX® PAD CONNECT is a sealing plate made of soft, skin-coloured medical silicone.

It is used to align the skin level of an irregular tracheostoma that is situated below the skin level.

The plano-convex shape of the LARYVOX® PAD CONNECT ensures high quality sealing.

Due to its material properties and ergonomic shape, the LARYVOX® PAD CONNECT adapts to the skin around the tracheostoma.

In addition to reduced leakage of secretion, the lateral flow of breathing air can also be minimised.

LARYVOX® PAD CONNECT has an integrated 22 mm connector so that a speaking valve/ HME filter cassettes can be fixated directly into the PAD.

LARYVOX® PAD CONNECT is particularly skin-friendly due to the material used.

No adhesive is required to fix the pad.

VI. INSTRUCTIONS

CAUTION!

When opening the packaging, make sure that the LARYVOX® PAD CONNECT is undamaged and that the indicated use-by date has not been exceeded.

Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or OPTIFAHL® stoma cleansing wipes (REF 33200).

Dry the skin.

The skin surrounding the tracheostoma must be clean, dry and free from grease.

After removing the protective foil backing ❶, place the LARYVOX® PAD CONNECT with the adhesive side centred over the tracheostoma ❷.

You can then insert a HME filter cassettes/speaking valve into the integrated 22 mm connector ❸/❹.

Before removing the LARYVOX® PAD CONNECT, the inserted aids must be removed.

To do this, fixate the LARYVOX® PAD CONNECT with two fingers and carefully remove the HME filter cassettes/ speaking valve from the pad ❺ to achieve a long retention time of the pad.

Remove the LARYVOX® PAD CONNECT by carefully removing the adhesive strips from the skin. ❻

CAUTION!

Do not use any fatty creams!

Otherwise, secure attachment of the product is no longer ensured.

VII. CLEANING

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. Please consult your physician.

Cleaning Steps

Rinse the LARYVOX® PAD CONNECT thoroughly under lukewarm running water.

Then allow LARYVOX® PAD CONNECT to dry in a clean environment.

VIII. STORAGE

Currently not used and cleaned LARYVOX® PAD CONNECT should be stored in a clean plastic container in a dry environment away from dust, sunlight and/or heat.

Still packaged replacement pads should be stored in a dry environment away from sunlight and/or heat.

We highly recommend keeping a replacement pad on hand to ensure continuous supply.

IX. PERIOD OF USE

After 29 days, the LARYVOX® PAD CONNECT should be disposed of properly and replaced with a new product.

Deviations from the expected wearing period occur as a result of individual wearing conditions and skin properties.

CAUTION!

Should the material show cracks or other damage, the product must be replaced immediately.

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Product specifications described herein are subject to change without notice. LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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